

HHS Public Access

Author manuscript *Nat Methods*. Author manuscript; available in PMC 2022 March 30.

Published in final edited form as:

Nat Methods. 2021 October ; 18(10): 1233-1238. doi:10.1038/s41592-021-01277-2.

Image-Guided MALDI Mass Spectrometry for High-Throughput Single-Organelle Characterization

Daniel C. Castro^{1,4}, Yuxuan Richard Xie^{4,5}, Stanislav S. Rubakhin^{2,3,4}, Elena V. Romanova^{2,3,4}, Jonathan V. Sweedler^{1,2,3,4,5,*}

¹Department of Molecular and Integrative Physiology, University of Illinois at Urbana-Champaign, Urbana, IL USA

²Department of Chemistry, University of Illinois at Urbana-Champaign, Urbana, IL USA

³Neuroscience Program, University of Illinois at Urbana-Champaign, Urbana, IL USA

⁴Beckman Institute for Advanced Science and Technology, University of Illinois at Urbana-Champaign, Urbana, IL USA

⁵Department of Bioengineering, University of Illinois at Urbana-Champaign, Urbana, IL USA

Abstract

Peptidergic dense-core vesicles are involved in packaging and releasing neuropeptides and peptide hormones – critical processes underlying brain, endocrine and exocrine function. Yet, the heterogeneity within these organelles, even for morphologically defined vesicle types, is not well characterized because of their small volumes. We present image-guided, high-throughput mass spectrometry-based protocols to chemically profile large populations of both dense-core vesicles and lucent vesicles for their lipid and peptide contents, allowing observation of the chemical heterogeneity within and between these two vesicle populations. The proteolytic processing products of four prohormones are observed within the dense-core vesicles, and the mass spectral features corresponding to the specific peptide products suggest three distinct dense-core vesicle populations. Significant differences in the lipid mass range are observed between the dense-core and lucent vesicles. These single-organelle mass spectrometry approaches are adaptable to characterize a range of subcellular structures.

Organelles are one of the smallest structural units that influence the functional, morphological, and biochemical characteristics of different cell types. Chemical analysis of

Author contributions

Competing interests

Additional information

Users may view, print, copy, and download text and data-mine the content in such documents, for the purposes of academic research, subject always to the full Conditions of use: https://www.springernature.com/gp/open-research/policies/accepted-manuscript-terms

^{*}Corresponding author, jsweedle@illinois.edu.

J.V.S. conceived the project idea. S.S.R. and D.C.C. designed and performed the sample preparation approach. D.C.C. acquired microscopy images, developed the image-processing approaches and acquired the high-throughput MALDI MS data. Y.R.X. performed data analysis and contributed to the figures. E.V.R. performed LC-MS/MS peptide sequencing and bioinformatics. All authors wrote and approved the final manuscript. J.V.S. acquired funding and provided direction throughout the project.

The authors declare no competing interests.

Supplementary information is available for this paper.

individual organelles is challenging due to their attoliter volumes, the wide dynamic range of analyte concentrations, and the need for sophisticated isolation procedures, thereby limiting our understanding of their chemical heterogeneity. Mass spectrometry imaging has begun being used for both cellular and some subcellular analyses in discovery-based studies^{1–3} but is limited in throughput and spatial resolution for organelle measurements. Alternatively, both cells^{4–6} and organelles⁷ of interest can be isolated and placed on a glass slide for subsequent mass spectrometry (MS) measurement. Though we used this approach to assay the peptides within individual organelles⁷, it was limited in throughput and difficult to automate due to the manual positioning of each individual organelle prior to measurement.

A recent enhancement to single-cell measurement⁸ involves scattering the cells of interest onto a microscope slide, determining their locations via fluorescence microscopy, and then targeting selected locations with matrix-assisted laser desorption/ionization (MALDI) MS, allowing tens of thousands of cells to be assayed. Here we adapted the single-cell approach to single organelles, which required three enhancements: (1) improved object targeting approaches, (2) optimized analyte detection using high-resolution mass spectrometry and (3) unsupervised data analysis workflows to characterize organelle heterogeneity. These advances allowed us to use MALDI Fourier-transform ion cyclotron resonance MS for the high-throughput simultaneous detection of both peptide and lipid species in 0.5 to 2 µm diameter dense-core vesicles (DCVs) and electron lucent vesicles (LVs) isolated from the exocrine atrial gland (AG) (Fig. 1a-f) and red hemiduct (RH) of Aplysia californica, respectively. By developing single-organelle sampling techniques that can be conducted in a high-throughput manner, we observed subtypes of DCVs defined by their overlapping but distinct peptide content. We also identified a peptide prohormone not previously known to localize within DCVs⁹⁻¹¹, as well as large differences in the contents between the DCVs and LVs. Our method revealed post-translational proteolytic processing of AGPB1 (XP 012945142.1), AGPA1 (XP 012945143.1) and AGPA2 (XP 012945134.1) prohormones, with peptide structures validated by liquid chromatography-tandem mass spectrometry (LC-MS/MS). Additionally, a novel prohormone, atrial gland Peptide D, was characterized (Supplementary Figs. 1 and 2; Supplementary Table 1). Importantly, our workflow is extendable to multiple imaging modalities such as scanning electron microscopy (Supplementary Fig. 3) and instrumentation that uses a piezo linear stage and camera, providing an avenue for the analysis of targets that are smaller than the wavelength of light.

Results

Image-guided MALDI MS approach.

In order to adapt single-cell approaches to objects that are 10-fold smaller in diameter than cells measured previously, we created a method to enhance identification of micronsized object locations while improving analyte detectability by minimizing the number of interfering chemicals that could reduce the ionization efficiency of our analytes of interest. Traditionally, single-cell targeting requires chemical labeling to mark objects of interest, but introduction of exogenous chemicals during single-organelle sample preparation is problematic due to the 1000-fold less material present. Paraformaldehyde fixation is a

hallmark step in most staining protocols, but paraformaldehyde fixation crosslinks DNA and peptides/proteins, rendering them poorly ionizable and undetectable by MS. By measuring the pixel properties produced by DCVs under brightfield conditions, we compiled sequences of image algorithms (Supplementary Fig. 4) that selectively identify the pixels produced by DCVs versus the background (Fig. 1b-d), allowing identification of the spatial locations of individual DCVs without the need for chemical labeling. The axial resolution of light microscopy limits this approach to the analysis of vesicles that are approximately 500 nm in diameter or larger. However, non-optical approaches such as scanning electron microscopy can be applied, providing a resolution down to 10 nm, allowing targeting of nanometer-sized objects (Supplementary Fig. 3). A three-step approach was developed for isolation of a representative number of vesicles without loss of chemical detail, which is critical for the high-throughput preparation of organelle samples. In this approach, a volatile, isosmotic-ammonium acetate buffer was deposited onto an indium tin oxide (ITO) coated glass slide, into which a small aliquot of artificial sea water solution containing the vesicles was deposited. The vesicles were allowed to sediment and adhere to the glass surface before the solution was aspirated, leaving a large number of visually intact vesicles distributed across the slide. After identifying the vesicles on the slide using the described pipeline, an image mask of the identified objects was created. The image mask marks the vesicles with a maximal pixel intensity and forces the remaining pixels to a zero-pixel intensity, creating a binary image output (Fig. 1d). With the vesicles marked by clusters of maximal pixel intensity values, the binary image created from the masked image output can be used for identification of the spatial locations of vesicles on the slide for subsequent analysis via MALDI MS. A 200 µm distance filter was applied to remove objects from the target list that are closer than 200 µm to each other, ensuring the 100 µm diameter MALDI laser spot size did not overlap with multiple vesicles (Supplementary Fig. 5).

High-throughput DCV detection and characterization.

A total of 598 DCVs were analyzed and after cross-refencing the single-DCV spectra with the LC-MS/MS data, over 50 mature full-length peptides were assigned by peptide mass fingerprinting to eight known prohormones and one novel prohormone (Supplementary Table 2). As a result, the peptides AGPA2 [36–69], AGPB1 [85–118], AGPA1 [117–152], AGPA1 [156–173] (or identical peptide AGPB1 [205–222]), N-terminal peptide AGPA1 [22–34], and the truncated N-terminal peptide AGPB1 [71–81] lacking the C-terminal Thr-Asn⁷ were detected, resulting in 39% and 28% coverage of AGPA1 and AGPB1 prohormones (Fig. 1g,h), respectively. Lastly, by using an *m/z* range of 150–4,500, we were able to perform simultaneous detection of both peptide and lipid species in individual DCVs. Three phosphatidylcholine (PC) lipid species detected by mass-match assignment include: PC(18:1/16:0), PC(18:1e/16:0) and PC(18:1e/18:1), with respective assignment mass errors of -1.62 ± 0.65 ppm, -2.59 ± 0.06 ppm and -1.55 ± 0.21 ppm (Fig. 2e). The three detected lipids were previously identified in *Aplysia californica* neurons using MALDI MS¹².

Differential packaging of peptides from the same prohormone into different vesicles was shown in the bag cell neurons of *Aplysia californica* using immunogold staining approaches¹³, although this has not been demonstrated in other *Aplysia californica* cell types. Here, we considered whether the different peptide complements may be

from differential proteolytic processing of AGPB1, AGPA1 and AGPA2 prohormones. Therefore, CX decomposition, an unsupervised statistical approach, was applied to select statistically important features while removing redundant and uninformative ones for improved downstream data interpretation (Supplementary Figs. 6 and 7). CX decomposition based on statistical leverage scores was used to select the top 200 features¹⁴ (Fig. 2a,b; Supplementary Fig. 8). Next, using the selected features, the DCVs were partitioned into three groups using the k-means clustering algorithm, for which the number of groups was determined by calculating the Within Cluster Sum of Squares (Fig. 2c). Significantly different features were obtained by performing a Wilcoxon rank-sum test for each individual cluster against the others (Supplementary Fig. 9). The violin plots (Fig. 2d) highlight selected mass spectral features and their distribution in the different DCV cluster types where known AG peptides are labeled. The N-terminal peptides, AGPB1 [71-81] and AGPA1 [22–34], and the novel Peptide D [63–75] (Fig. 2f), are shown to be defining features in Cluster 1 when compared to Cluster 2, whereas Cluster 3 demonstrates localization of multiple expected AG peptides from the prohormone C-termini, including AGPA1 [156–173], Peptide D [63–75], AGPB1[85–118] and a peptide cluster belonging to Peptide D [132-162] (Supplementary Fig. 10). The data indicate that the N- and C-terminal peptides from these two prohormones are differentially packaged into a subset of DCVs, whereas the bioactive peptides that initiate the egg-laying process¹⁵, AGPA2 [36–69] and AGPB1 [85-118], are predominantly packaged into a separate subset of vesicles. An early pioneering study that quantified the amounts of AGPA1 [22-34], AGPA2 [36-69] and AGPB1 [85–118] in isolated intact granules from AG homogenates also suggested the possibility of packaging of AGPA1 [22–34] into a distinct class of secretory vesicle¹⁵. Perhaps this differential packaging of N-terminal products in the AG would be akin to what is also observed in the bag cell neurons where N- and C-terminal peptides are differentially packaged into two discrete vesicle classes that are targeted to different cellular locations¹³.

Mass spectrometric characterization of lucent vesicles.

To demonstrate the suitability of our method to assay a variety of organelles, we extended it to another morphologically distinct vesicle type, LVs, from the adjacent RH of Aplysia californica. Simple modification to the image-processing pipeline allows direct targeting of the LVs based on their different morphological characteristics (Fig. 3a). The LVs of the RH have a similar diameter of 0.5 to 2 µm but vary dramatically in ultrastructure; the LVs contain a double membrane structure, with the inner membrane forming cristae-like structures that are analogous to those found within the mitochondrion^{9,16}. In contrast, DCVs from the AG contain a single membrane and exhibit a solid dense-core compared to the more relatively fluid interior of the LVs from the RH^{9,16}. Due to the nonpeptidergic content of the LVs and the clear differences in ultrastructure between the LVs and DCVs, a machine learning model was trained for the differentiation between LVs and DCVs, focusing on the m/z range of 500–1,100. t-distributed stochastic neighborhood embedding (t-SNE) was initially applied using the mass spectral features to visualize the data in a low-dimensional space (Figure. 3b). Next, gradient boosting trees were trained with a 3-fold validation, obtaining an accuracy of $98.6 \pm 0.78\%$ for the classification between LVs and DCVs (Fig. 3c). The most important features for the classification between LVs and DCVs (Fig. 3d) were selected via Shapley additive explanations (SHAP) through a previously described

method ⁶, where a total of 97 features with non-zero mean SHAP values from the model output were selected, with 36 features putatively identified by mass-match assignment (Fig. 3e). The model provided an improved classification performance through retraining on the selected features (Supplementary Fig. 11), displaying the differential lipid profiles between the two vesicle classes (Supplementary Figs. 12, 13 and 14). The data show the strong presence of multiple sterol lipid species to be of importance to their classification. The absence of sterol lipid species correlates with uniformly curved membranes^{17,18}. Whereas when present in sufficient concentration, sterol lipid species will distribute asymmetrically in the lipid bilayer causing large, flat membrane regions that are separated by sharp curves in the membrane, creating negative membrane curvature. Interestingly, negative membrane curvature would promote the formation of the cristae-like invaginations that can be observed in the LVs^{16,18}.

Discussion

As the spatial resolution of MS advances from the single-cell level to include the singleorganelle domain, new analytical approaches are being developed to unravel the complex microenvironments housed within single cells^{19,20,21}. Capillary micro-sampling techniques coupled to electrospray ionization (ESI) MS have been developed for the identification of peptides in individual cells and in some subcellular compartments of cells^{11,22}. Though capillary micro-sampling ESI-MS techniques can identify peptides in subcellular compartments^{23,24}, they require fine-tuning of the capillary position and controlling the volume extracted. These requirements are critical to single-organelle extraction and pose significant challenges to repeated sampling and throughput due to potential contamination from the cellular membrane, cytoplasm and other subcellular organelles²². Alternatively, our image-guided MALDI MS approach provides direct single-organelle targeting using point-based image registration.

In conventional mass spectrometry imaging, as spatial resolution increases, so does acquisition time. Therefore, to maintain the throughput required for characterizing large populations of individual organelles, we are not imaging slides containing the vesicles. Instead, the organelles themselves are isolated and selectively targeted using the MALDI laser. Our high-throughput, label-free, single-organelle approach enables simultaneous characterization of both the peptide and lipid content of hundreds of intact 0.5 to 2 μ m diameter DCVs and LVs. Because our approach uses a combination of freeware and instrumentation that is commonly available to researchers, it is adaptable to many experimental designs, including other organelle types such as microvesicles, which can range from 100 nm up to 1 μ m in diameter²⁵.

Online Methods

General.

All chemicals were obtained from Sigma-Aldrich (St. Louis, MO) unless specified otherwise.

Vesicle isolation.

Aplysia californica (150–200 g) were obtained from National Resource for Aplysia (Miami, FL) and kept in 14 °C aquarium filled with Instant Ocean (Aquarium Systems Inc., Mentor, OH). Animal euthanasia was performed in accordance with the AVMA Guidelines for the Euthanasia of Animals: 2020 Edition (Section S6.3.1.1 Non-inhaled agents for immersion). Animals were anesthetized using a 333 mM MgCl₂ solution injected into the body cavity (50% volume/body).

For DCV isolation, the AG was isolated by manual dissection, placed into a microcentrifuge tube containing 1 mL of artificial sea water (ASW). Due to the holocrine release mechanism, DCV secretion from the AG was induced by gentle trituration with a polypropylene Pasteur pipette, releasing intact DCVs into the ASW solution. 100 μ L of 500 mM ammonium acetate solution was pre-spotted onto an ITO-unpolished float glass slide, $R_s = 70-100 \Omega$ (Delta Technologies, Loveland, CO). 50 μ L of vesicle-containing ASW solution was then pipetted into the 100 μ L of ammonium acetate solution on the ITO-glass slide for a total volume of 150 μ L. The resulting vesicle solution, leaving DCVs seeded across the slide. All steps were visually monitored using an inverted microscope. Each slide held DCVs from three animals (biological replicates) and a total of three slides (technical replicates) were prepared. For LV isolation, the procedure used for the DCVs was repeated using red hemiduct tissue, with the other steps being the same as for the DCV isolation.

Vesicle imaging.

Brightfield images were acquired on an Axio Imager M2 (Zeiss, Jena, Germany) equipped with an AxioCam ICc 5, a $.63 \times$ camera adaptor and TL VIS-LED lamp. Images were acquired in mosaic mode using a $10 \times$ objective with 30% tile overlap. The resulting tiles were stitched before exporting in TIFF-file format using ZEN 2.0 Pro edition (Zeiss, Jena, Germany) software.

Scanning electron microscopy images were acquired using an FEI Quanta FEG 450 environmental scanning electron microscope (FEI, Hillsboro, OR). Images were acquired using an accelerating voltage of 10 kV, dwell time of 10 µs and working distance of 6.6 mm.

Image-processing for vesicle recognition in microMS.

We utilized biological-image analysis software, including ImageJ²⁶ and CellProfiler²⁷, for vesicle pixel recognition on the brightfield microscopy image. Importantly, other imageprocessing software can be used for this process as well. CellProfiler pipeline modules were created to selectively mask the DCVs and LVs on the brightfield microscopy image to remove the background and retain just information on the vesicles of interest as the foreground. Masking the identified vesicles creates a binary-image marking the pixel locations of vesicles on the glass slide. The microMS software was used to translate the pixel locations on microscopy images to physical coordinates of the mass spectrometer's stage. A 200 μ m distance filter was then applied in microMS (removing vesicles located closer than 200 μ m to each other from target list). Alternatively, ImageJ was used as a

simple thresholding strategy by measuring the difference between the vesicle (foreground) and background pixel intensity, which allows a threshold to be set, leaving only objects of interest for recognition in microMS.

Matrix application.

Matrix deposition for MALDI MS analysis was performed using a glass sublimation apparatus (Wilmad-LabGlass, Vineland, NJ) filled with 2,5-dihydroxybenzoic acid as the MALDI matrix. The slide was attached to the cold-finger and vacuum was created using a rotary vane pump (Edwards Vacuum, model E2M30, Soquel, CA). The sublimation apparatus was placed on a sand bath preheated to 150 °C for 8 min. Matrix deposition was followed by recrystallization using 5% methanol. Recrystallization was performed using a 100 \times 15 mm polystyrene Petri dish as a recrystallization chamber. The ITO-glass slide was attached to the top of the recrystallization chamber and a filter paper (Whatman® Grade 1 Qualitative Filter Paper, Thermo Fisher Scientific, Waltham, MA) was wetted with 1 mL of 5% methanol. The chamber was sealed using tape and placed in an oven at 85 °C for 1.5 min. After removal of the recrystallization chamber from the oven, the slide was immediately removed from the recrystallization chamber and allowed to dry in a nitrogen chamber until analysis.

MALDI MS measurements.

High-throughput single-DCV and single-LV analysis were performed on a SolariX XR 7T Fourier-transform ion cyclotron resonance mass spectrometer equipped with an APOLLO II dual MALDI/ESI source (Bruker Corp., Billerica, MA) using an m/z range of 150–4,500. Data were acquired at 1 M giving a mass resolution of 107,000 at m/z 535 and 19,070 at m/z 3,922, yielding a transient length of .721s. The instrument was operating in positive-mode using a Smartbeam-II UV laser (Bruker Corp.) set to 'Ultra mode', which yields a 100 μ m diameter laser footprint. Each MALDI acquisition consisted of two accumulations comprised of 400 laser shots each, at a frequency of 1,000 Hz. DCV and LV stage coordinates and geometry files were generated using microMS as previously described⁸.

Data preprocessing.

Data preprocessing was performed using Compass Data Analysis 4.4.2 (Bruker Corp.) and MATLAB 2018b (MathWorks, Inc., Natick, MA). Internal calibration using the exact mass of AG peptides was performed using AGPB1 [71–81] (m/z 1,221.6878), AGPA1 [22–34] (m/z 1,396.7225), AGPA1 [156–173] (m/z 1,908.8761), AGPA2 [33–69] (m/z 3,922.9478) and AGPB1 [85–118] (m/z 4,031.0053). Peak picking and peak export for statistical analysis in MATLAB 2018b was set to S/N 5 with a relative intensity threshold of 0.01%. A non-uniform bin width was used for mass spectral alignment. For DCV data analysis, mass features were truncated at m/z 1,100 for downstream data analysis of only AG peptides. For LV and DCV vesicle classification, mass features in m/z 500–1,100 were selected for downstream data analysis. Internal calibration using AG peptides was not performed for the LV and DCV classification tasks. Using the target list provided by microMS, the pixel coordinates were used to find and crop the locations of individual vesicles across the original microscopy image. The corresponding mass spectra were matched with the appropriate vesicle for visual evaluation of corresponding single vesicles.

CX decomposition and statistical analysis.

Unsupervised approaches such as principal component analysis (PCA) and more generally, matrix decomposition or factorization, are valuable data analysis tools to enable field-specific interpretation of high-dimensional datasets. However, the interpretation is limited due to the complicated eigenspace obtained from PCA, deterring our further understanding of the feature space of the data matrix. CX matrix decomposition is designed to obtain a low-rank approximation of the data in terms of actual rows or columns¹⁴. Given an $m \times n$ data matrix A, the algorithm decomposes it into an $m \times c$ matrix C and a $c \times n$ matrix X, where C is expressed by c number of column vectors of the original data. The statistical leverage scores are used to rank and select the columns of C from A, which can be obtained by

$$l_j = \sum_{i=1}^k v_{ji}^2$$

where I_j is the leverage score for the *j*th column/feature, v_i is the right singular vectors obtained by the singular value decomposition, and *k* is the rank to be selected. *X* is then determined by minimizing the error:

$$\min_{X} \|A - CX\|_F$$

The reconstruction error evaluation and the rank k selection are provided as Supplementary Fig. 8. Based on the evaluation, the top 200 features are selected to form the columns of the matrix *C*. K-means clustering, and the Wilcoxon rank-sum test, were performed using the Python-based open source package SCANPY²⁸. The stacked violin plots, shown in Fig. 2d, of the normalized signal intensities represent the identified peak features across the 3 clusters, with the vesicle type (or cluster) assignments obtained by K-means clustering. The y-axis is the root-mean-squared normalized peak intensity and each violin contains the distribution of the normalized intensities of the corresponding features in the x-axis.

Peptide sequencing by LC-MS/MS.

Peptide extracts (n = 3) were obtained by manually grinding entire AG tissue in 500 µL of acidified methanol followed by evaporation and reconstitution of each extract in 0.1% formic acid. A nanoElute (Bruker Corp.) ultra-high-pressure nano-flow chromatography system was coupled to a trapped ion mobility - quadrupole time-of-flight mass spectrometer (timsTOF Pro, Bruker Corp.) with a CaptiveSpray nano-electrospray ion source (Bruker Corp.) equipped with an external column oven. Mobile phases A and B were water with 0.1% formic acid (v/v) and acetonitrile with 0.1% formic acid (v/v), respectively. Samples were loaded onto a pre-column peptide trap (AcclaimTM PepMapTM 100 C18, 1×5 mm, 5 µm particle size, Thermo Fisher Scientific) using solvent A for off-line desalting. Next, the trap was placed in-line with the analytical column and peptide separation was performed at 40 °C with a uniform flow of 300 nL/min on a C18 ReproSil AQ column (Bruker FIFTEEN, P/N#1842621: 150 mm × 75 µm, 1.9 µm particle size, pore size 120A) equilibrated at 2%

B. A linear gradient of solvent B was applied as follows: 2-10% within 5 min, 10-50% in the next 115 min, followed by a washing step at 95% B and re-equilibration, during which data collection was not performed. The mass spectrometer was operated in parallel accumulation-serial fragmentation (PASEF) mode for peptide sequencing. The mass range for precursor ion was set to m/z 100–1,700, ion mobility 1/k0 range 0.6–1.6 Vs/cm². Fragmentation was performed with 10 PASEF scans, cycle time of 1.1 sec, during which collision energy varied linearly between 20 and 59 eV depending on precursor 1/k0 value within the set range. Active dynamic exclusion of precursor ions was set to 0.4 min.

Peptide identification by bioinformatics.

MS raw files were processed with PEAKS Online²⁹ (Bioinformatics Solution Inc., Waterloo, Canada) using the DeNovo, DB and PTM protocols, sequentially. Peptide sequence tags obtained by the DeNovo process (80% ALC score cut off) were searched against the *Aplysia* RefSeq database available from the NCBI (GCF_000002075.1). The protein database was filtered to include proteins up to 1,000 amino acids long. Search parameters for DB included: parent mass error tolerance 20.0 ppm, fragment mass error tolerance 0.03 Da, no enzyme, digest mode unspecific. Next, PTMs were identified by searching the data for amidation, acetylation (K) and acetylation N-terminus, pyro-glutamylation from (Q) and (E), phosphorylation (STY) and half-disulfide bridge. The false discovery rate was determined by decoy fusion method and the threshold set to 1% for peptides. Proteins with –10lgP>20 and at least one unique peptide are reported.

Vesicle classification through machine learning.

We adapted a previously described machine learning strategy⁶ for vesicle type classification to differentiate between DCVs and LVs based on their lipid contents. Features in the m/z range of 500–1,100 were used for the classification task. Gradient boosting trees were trained with a 3-fold validation; in each fold of the model performance, the metrics were computed to obtain classification accuracies, confusion matrices and receiver operating characteristic (ROC) curves (Supplementary Fig. 11). The most contributing features to the classification task were selected via SHAP, a game theory approach for model explanations, which were obtained through a Python implementation of SHAP. A total of 97 features with non-zero mean SHAP values from the output of trained models were selected and 36 were putatively annotated by searching against the LIPID MAPS³⁰ database with a 7-ppm tolerance. Models were then retrained with the SHAP-selected features as well as the annotated features to verify the discriminative ability of those features. t-SNE using the cosine distance was used to visualize the lipid differences between DCVs and LVs in a low-dimensional space.

Reporting Summary

Further information on research design is available in the Nature Research Reporting Summary linked to this article.

Data availability

The data that supports the findings of this study are publicly available via the Illinois Data Bank (https://doi.org/10.13012/B2IDB-5949772_V1).

Code availability

The code used in this study is available on GitHub (https://github.com/richardxie1119/ DCV DA).

Supplementary Material

Refer to Web version on PubMed Central for supplementary material.

Acknowledgments

This work was supported by funding from the National Institute on Drug Abuse under Award No. P30DA018310 (to J.V.S.) and the National Human Genome Research Institute under award No. R01HG010023 (to J.V.S.). The content is solely the responsibility of the authors and does not necessarily represent the official views of the awarding agencies.

References

- Lovric J et al. Nano Secondary Ion Mass Spectrometry Imaging of Dopamine Distribution Across Nanometer Vesicles. ACS Nano 11, 3446–3455, (2017). (doi:10.1021/acsnano.6b07233) [PubMed: 27997789]
- Ibanez AJ et al. Mass spectrometry-based metabolomics of single yeast cells. Proc. Natl. Acad. Sci. U. S. A 110, 8790–8794, (2013). (doi:10.1073/pnas.1209302110) [PubMed: 23671112]
- Kompauer M, Heiles S & Spengler B Atmospheric pressure MALDI mass spectrometry imaging of tissues and cells at 1.4-mum lateral resolution. Nat. Methods 14, 90–96, (2017). (doi:10.1038/ nmeth.4071) [PubMed: 27842060]
- Neumann EK, Ellis JF, Triplett AE, Rubakhin SS & Sweedler JV Lipid Analysis of 30000 Individual Rodent Cerebellar Cells Using High-Resolution Mass Spectrometry. Anal. Chem 91, 7871–7878, (2019). (doi:10.1021/acs.analchem.9b01689) [PubMed: 31122012]
- Neumann EK, Do TD, Comi TJ & Sweedler JV Exploring the Fundamental Structures of Life: Non-Targeted, Chemical Analysis of Single Cells and Subcellular Structures. Angew. Chem. Int. Ed. Engl 58, 9348–9364, (2019). (doi:10.1002/anie.201811951) [PubMed: 30500998]
- Xie YR, Castro DC, Bell SE, Rubakhin SS & Sweedler JV Single-Cell Classification Using Mass Spectrometry through Interpretable Machine Learning. Anal. Chem 92, 9338–9347, (2020). (doi:10.1021/acs.analchem.0c01660) [PubMed: 32519839]
- Rubakhin SS, Garden RW, Fuller RR & Sweedler JV Measuring the peptides in individual organelles with mass spectrometry. Nat. Biotechnol 18, 172–175, (2000). (doi:10.1038/72622) [PubMed: 10657123]
- Comi TJ, Neumann EK, Do TD & Sweedler JV microMS: A Python Platform for Image-Guided Mass Spectrometry Profiling. J. Am. Soc. Mass Spectrom 28, 1919–1928, (2017). (doi:10.1007/ s13361-017-1704-1) [PubMed: 28593377]
- Ciobanu L et al. Characterization of the physicochemical parameters of dense core atrial gland and lucent red hemiduct vesicles in Aplysia californica. Anal. Chem 76, 2331–2335, (2004). (doi:10.1021/ac035346h) [PubMed: 15080745]
- Nagle GT, Painter SD, Blankenship JE & Kurosky A Proteolytic processing of egg-laying hormone-related precursors in Aplysia. Identification of peptide regions critical for biological activity. J. Biol. Chem 263, 9223–9237, (1988). [PubMed: 3379066]

- Lillard SJ et al. Separation and characterization of amines from individual atrial gland vesicles of Aplysia californica. Anal. Chem 70, 3517–3524, (1998). (doi:10.1021/ac980462j) [PubMed: 9737203]
- Passarelli MK, Ewing AG & Winograd N Single-cell lipidomics: characterizing and imaging lipids on the surface of individual Aplysia californica neurons with cluster secondary ion mass spectrometry. Anal. Chem 85, 2231–2238, (2013). (doi:10.1021/ac303038j) [PubMed: 23323749]
- Conn PJ & Kaczmarek LK The bag cell neurons of Aplysia. A model for the study of the molecular mechanisms involved in the control of prolonged animal behaviors. Mol. Neurobiol 3, 237–273, (1989). (doi:10.1007/BF02740607) [PubMed: 2698177]
- Yang J, Rubel O, Prabhat, Mahoney, M. W. & Bowen, B. P. Identifying important ions and positions in mass spectrometry imaging data using CUR matrix decompositions. Anal. Chem 87, 4658–4666, (2015). (doi:10.1021/ac5040264) [PubMed: 25825055]
- Nagle GT, Van Heumen WR, el-Hamzawy MA & Kurosky A Purification and characterization of Aplysia atrial gland secretory granules containing egg-laying prohormone-related peptides. Peptides 15, 101–108, (1994). (doi:10.1016/0196-9781(94)90177-5) [PubMed: 8015966]
- Beard M, Millecchia L, Masuoka C & Arch S Ultrastructure of secretion in the atrial gland of a mollusc (Aplysia). Tissue Cell 14, 297–308, (1982). (doi:10.1016/0040-8166(82)90028-3) [PubMed: 7112536]
- 17. Barenholz Y Cholesterol and other membrane active sterols: from membrane evolution to "rafts". Prog. Lipid Res 41, 1–5, (2002). (doi:10.1016/s0163-7827(01)00016-9) [PubMed: 11694266]
- Yesylevskyy SO, Demchenko AP, Kraszewski S & Ramseyer C Cholesterol induces uneven curvature of asymmetric lipid bilayers. ScientificWorldJournal 2013, 965230, (2013). (doi:10.1155/2013/965230) [PubMed: 23766730]
- Mathew DG et al. Electrochemical Detection of Tumor-Derived Extracellular Vesicles on Nanointerdigitated Electrodes. Nano Lett 20, 820–828, (2020). (doi:10.1021/acs.nanolett.9b02741) [PubMed: 31536360]
- Zhou S et al. Integrated Microfluidic Device for Accurate Extracellular Vesicle Quantification and Protein Markers Analysis Directly from Human Whole Blood. Anal. Chem 92, 1574–1581, (2020). (doi:10.1021/acs.analchem.9b04852) [PubMed: 31779307]
- Théry C, Amigorena S, Raposo G & Clayton A Isolation and characterization of exosomes from cell culture supernatants and biological fluids. Curr. Protoc. Cell Biol Chapter 3, Unit 3 22, (2006). (doi:10.1002/0471143030.cb0322s30)
- Zhang L et al. Subcellular Peptide Localization in Single Identified Neurons by Capillary Microsampling Mass Spectrometry. Sci. Rep 8, 12227, (2018). (doi:10.1038/s41598-018-29704-z) [PubMed: 30111831]
- Zhang L & Vertes A Single-Cell Mass Spectrometry Approaches to Explore Cellular Heterogeneity. Angew. Chem. Int. Ed. Engl 57, 4466–4477, (2018). (doi:10.1002/anie.201709719) [PubMed: 29218763]
- Hu K, Nguyen TDK, Rabasco S, Oomen PE & Ewing AG Chemical Analysis of Single Cells and Organelles. Anal. Chem 93, 41–71, (2021). (doi:10.1021/acs.analchem.0c04361) [PubMed: 33284577]
- Stahl AL, Johansson K, Mossberg M, Kahn R & Karpman D Exosomes and microvesicles in normal physiology, pathophysiology, and renal diseases. Pediatr. Nephrol 34, 11–30, (2019). (doi:10.1007/s00467-017-3816-z) [PubMed: 29181712]

References

- Rueden CT et al. ImageJ2: ImageJ for the next generation of scientific image data. BMC Bioinformatics 18, 529, (2017). (doi:10.1186/s12859-017-1934-z) [PubMed: 29187165]
- Jones TR et al. CellProfiler Analyst: data exploration and analysis software for complex imagebased screens. BMC Bioinformatics 9, 482, (2008). (doi:10.1186/1471-2105-9-482) [PubMed: 19014601]
- Wolf FA, Angerer P & Theis FJ SCANPY: large-scale single-cell gene expression data analysis. Genome Biol 19, 15, (2018). (doi:10.1186/s13059-017-1382-0) [PubMed: 29409532]

- 29. Ma B et al. PEAKS: powerful software for peptide de novo sequencing by tandem mass spectrometry. Rapid Commun. Mass Spectrom 17, 2337–2342, (2003). (doi:10.1002/rcm.1196) [PubMed: 14558135]
- Fahy E et al. Update of the LIPID MAPS comprehensive classification system for lipids. J. Lipid Res 50 Suppl, S9–14, (2009). (doi:10.1194/jlr.R800095-JLR200) [PubMed: 19098281]

Castro et al.

Fig. 1:

High-throughput workflow for label-free single-DCV targeting and MS analysis. (a) Schematic of MALDI MS workflow for high-throughput single-DCV measurements. (b) Brightfield image of DCVs distributed on a glass slide. (c) Identification of primary objects (DCVs) using image-processing software. (d) Masked image output of identified DCVs. The colored spots represent primary objects (DCVs) recognized as 'foreground' and are marked with a maximal pixel intensity value. Anything not identified as an object is treated as 'background' and is set to a zero-pixel intensity value. (e) Brightfield image

of mechanically induced DCV release from the AG. (f) DCVs plated on a glass slide for relative DCV density estimation using brightfield microscopy. Each slide held DCVs from three animals (biological replicates) and a total of three slides (technical replicates) were prepared, where 598 DCVs were measured. (g) Mass spectra demonstrating the coverage of AG peptides detected in single-DCV measurements. (h) AGPA1 (XP_012945143.1) and AGPB1 (XP_012945142.1) prohormone sequences with corresponding MALDI MS-detected peptides *italicized* and font colored to match the annotated spectra in (g). AG peptide assignments were validated using LC-MS/MS and performed on AG extracts (n = 3). The peptides detected by LC-MS/MS are underlined in black.

Castro et al.

Page 15

Fig. 2:

Chemical heterogeneity of DCV populations: unsupervised data analysis, simultaneous acquisition of analyte classes and detection of novel Peptide D. (a) Deterministic CX was used to select 200 mass spectral features to improve data interpretation. CX decomposition was applied for feature selection where the best low-rank approximation was provided through singular value decomposition with a rank parameter k = 150, which was determined by the reconstruction errors with different rank parameters. Rank k = 150 has a reconstruction error less than 25% of the original dataset. (b) 'Statistical leverage scores'

for all spectral features were computed and plotted against the m/z axis. Multiple detected AG peptides can be seen with high leverage scores and are annotated with blue triangles. The 'Statistical leverage score' for each feature indicates its influence on the best low-rank fit of the data matrix. Selecting the top 200 features ensured that CX decomposition closely approximated the best low-rank fit of the original data matrix while removing uninformative information. (c) Results of k-means clustering of the dataset containing 200 selected features with the highest statistical leverage score. (d) Violin plots visualizing a subset of the selected features and their relative distribution in each cluster, with known AG peptides marked. (e) Simultaneous detection of PC(18:1/16:0), PC(18:1e/16:0) and PC(18:1e/18:1) with their corresponding sodiated adducts annotated. (f) The novel prohormone, Peptide D, was identified by MALDI MS mass-match assignment of Peptide D[63–75] and Peptide D[132–162].

Castro et al.

Fig. 3:

Comparison between two distinct vesicle types, LVs and DCVs. (a, top) Digitally inverted brightfield image showing LVs distributed across the glass slide. (a, bottom) Modification of the image-processing pipeline allows multiple morphologically distinct objects to be identified and added or removed for subsequent analysis. Green areas represent objects accepted for analysis. Magenta outline represents an object (here a spermatozoon) removed from analysis. A small number of spermatozoa (annotated) typically collected during RH LV sample preparation were used to demonstrate the effectiveness of our object filtering approach. A total of 123 LVs were measured from the same three biological replicates used in the DCV isolation. b) t-SNE was performed on the initial dataset for visualization of LVs and DCVs using all mass spectral features. c) Confusion matrix of the prediction on the test data using a 3-fold validation. d) Representative LV mass spectrum (red) overlaid on a representative DCV mass spectrum (grey). Representative spectra were not preprocessed. A subset of mass spectral features determined important via SHAP are annotated in the LV mass spectrum. e) Mass spectral feature contribution plot showing the importance of features to the vesicle classification task across m/z 500–1,100 on the x-axis. The y-axis shows the mean absolute SHAP value for the corresponding m/z feature (normalized between 0 and 1). Grey dots represent features with elevated mean signal intensities in DCVs and pink dots represent features with elevated mean signal intensities in LVs. The different sized dots

represent the respective p-values (two-sided Wilcoxon rank-sum test) for each corresponding feature.