



HHS Public Access

Author manuscript

Bio Protoc. Author manuscript; available in PMC 2017 August 08.

Published in final edited form as:

Bio Protoc. 2017 June 20; 7(12): . doi:10.21769/BioProtoc.2335.

Optogenetic Stimulation and Recording of Primary Cultured Neurons with Spatiotemporal Control

J r mie Barral* and Alex D Reyes

Center for Neural Science, New York University, New York, USA

Abstract

We studied a network of cortical neurons in culture and developed an innovative optical device to stimulate optogenetically a large neuronal population with both spatial and temporal precision. We first describe how to culture primary neurons expressing channelrhodopsin. We then detail the optogenetic setup based on the workings of a fast Digital Light Processing (DLP) projector. The setup is able to stimulate tens to hundreds neurons with independent trains of light pulses that evoked action potentials with high temporal resolution. During photostimulation, network activity was monitored using patch-clamp recordings of up to 4 neurons. The experiment is ideally suited to study recurrent network dynamics or biological processes such as plasticity or homeostasis in a network of neurons when a sub-population is activated by distinct stimuli whose characteristics (correlation, rate, and, size) were finely controlled.

Keywords

Primary culture of neurons; Optogenetics; Patterned optical stimulation

Background

Optogenetics provide a mean to control neuronal activity with millisecond precision. However, neurons are often activated simultaneously either by flashes of light that activate the whole population synchronously or by a light whose intensity is temporally modulated over the whole field of view (Boyden *et al.*, 2005). Yet, several methods exist to modulate the light spatially and have been used to uncage glutamate (Nawrot *et al.*, 2009) or activate channelrhodopsin (ChR2) expressing neurons (Guo *et al.*, 2009) (for review of available methods to stimulate neurons with both spatial and temporal resolution see Anselmi *et al.*, 2015).

To gain spatial control of the stimulation, a first possibility is to use a laser and move its beam quickly over different positions. For example, uncaging glutamate at different dendritic locations has been achieved by deflecting a laser beam with acousto-optic deflectors (Shoham *et al.*, 2005). This strategy is likely viable only if we modulate the light intensity sufficiently slowly over a limited area. Alternatively, a spatial pattern of light can be achieved using phase or intensity light modulators. Holographic technique based on

*For correspondence: barral@cns.nyu.edu.

phase modulation permits to obtain an image in three dimensions with a good spatial precision but patterns can be displayed at a rate of only 100 Hz (Papagiakoumou *et al.*, 2010). If a two dimensional pattern is sufficient, intensity modulation can simply be obtained by placing a projector or an array of LEDs in the conjugated plane of the sample (Farah *et al.*, 2007; Guo *et al.*, 2009). This technique has the advantages of being easy to implement, can target many regions of interest simultaneously and has the fastest temporal resolution.

Here we took advantage of a fast video projector based on the workings of a Digital Micromirror Device (DMD). A LED light source is split by an array of micromirrors that can be controlled with sub millisecond precision in order to display any arbitrary pattern of light (Barral and Reyes, 2016). An image of the projector is focalized to the sample plane via a pair of lenses and the microscope objective. The DMD technology offers an unprecedented temporal precision that enables to display patterns at 1.44 kHz and even faster DMDs are now available. In our settings, the resulting pixel size ($2.2 \times 1.1 \mu\text{m}$) was sufficiently small to stimulate single neurons.

To activate a single neuron, we selected a region of interest of $\sim 30 \times 30 \mu\text{m}$, centered at the soma of the neuron of interest and sent a 5 msec pulse of light. By designing patterns that are projected onto the sample, we could target independently and simultaneously a large number of neurons (10 to 100 neurons). Stimulated neurons were both excitatory and inhibitory (expression of ChR2 under the Synapsin promoter) and were activated by Poisson spike trains. The rate and correlation of the spike stimuli were controlled by the experimenter (see Barral and Reyes, 2016). By recording from neurons that expressed ChR2, we verified that stimulated neurons responded faithfully to the light pulses. We then recorded concurrently the membrane potentials of up to 4 neurons in cell-attached and in whole-cell configurations to isolate the spiking activity and the postsynaptic inputs, respectively.

Materials and Reagents

A. For the neuronal culture

1. Round coverslips, 25 mm diameter, German glass (Electron Microscopy Sciences, catalog number: 72196-25)
2. Siliconized low-retention microcentrifuge tubes (1.5 ml) (Fisher Scientific, catalog number: 02-681-331)
3. Low-retention pipet tips 200 μl (Fisher Scientific, catalog number: 02-717-165)
4. Low-retention pipet tips 1,000 μl (Fisher Scientific, catalog number: 02-717-166)
5. Disposable Petri dishes (35 \times 10 mm) (Corning, Falcon[®], catalog number: 351008)

6. Stericup-GP sterile vacuum filtration system, 0.22 μm , polyethersulfone (EMD Millipore, catalog number: SCGPU05RE)
7. Syringe filters; MCE membrane; pore size: 0.22 μm (EMD Millipore, catalog number: SLGS033SS)
8. Sterile transfer pipets (Fisher Scientific, catalog number: 13-711-20)
9. Conical sterile polypropylene centrifuge tubes (Thermo Fisher Scientific, Thermo Scientific™, catalog number: 339650)
10. Postnatal (P0–P1) mice
11. Hydrochloric acid (HCl) (Sigma-Aldrich, catalog number: 320331)
12. Nitric acid (Sigma-Aldrich, catalog number: 258121)
13. 70% ethanol
14. Poly-L-lysine hydrobromide–mol wt 70,000–150,000 (Sigma-Aldrich, catalog number: P1274)
15. Sodium tetraborate decahydrate 99.5% (Sigma-Aldrich, catalog number: B9876-500G)
16. Boric acid (cell culture tested) (Sigma-Aldrich, catalog number: B9645-500G)
17. HBSS (10 \times), no calcium, no magnesium, no phenol red (Thermo Fisher Scientific, Gibco™, catalog number: 14185052)
18. Penicillin-streptomycin (10,000 U/ml) (Thermo Fisher Scientific, Gibco™, catalog number: 15140122)
19. HEPES (Thermo Fisher Scientific, Gibco™, catalog number: 15630080)
20. Agar (Sigma-Aldrich, catalog number: A1296-100G)
21. Sodium bicarbonate (NaHCO₃) (Fisher Scientific, catalog number: S233)
22. Sodium pyruvate solution (Sigma-Aldrich, catalog number: S8636-100ML)
23. Sodium hydroxide (NaOH) (Fisher Scientific, catalog number: S318)
24. L-cysteine hydrochloride monohydrate (Sigma-Aldrich, catalog number: C7880-500MG)
25. DNase I grade II, from bovine pancreas–100 mg (Roche Diagnostics, catalog number: 10104159001)
26. Papain from *Carica papaya*–10 ml 100 mg (Roche Diagnostics, catalog number: 10108014001)
27. Calcium chloride (CaCl₂) (Fisher Scientific, catalog number: C79)

28. Magnesium chloride ($MgCl_2$) (Sigma-Aldrich, catalog number: M8266)
 29. Trypsin inhibitor from chicken egg white (Sigma-Aldrich, catalog number: T9253-500MG)
 30. Albumin from bovine serum (powder, suitable for cell culture, = 96%) (Sigma-Aldrich, catalog number: A9418-5G)
 31. Trypan blue solution (0.4%, liquid, sterile-filtered, suitable for cell culture) (Sigma-Aldrich, catalog number: T8154-20ML)
 32. Neurobasal medium (1×) (Thermo Fisher Scientific, Gibco™, catalog number: 21103049)
 33. B-27 serum-free supplement containing vitamin A (50×), liquid (Thermo Fisher Scientific, Gibco™, catalog number: 17504044)
 34. GlutaMAX™ supplement (Thermo Fisher Scientific, Gibco™, catalog number: 35050061)
 35. AAV virus for ChR2 expression (University of North Carolina Vector Core Services, AAV2-hSyn-hChR2(H134R)-mCherry)
 36. Sodium chloride (NaCl) (Fisher Scientific, catalog number: BP358)
 37. D-glucose (Fisher Scientific, catalog number: BP350)
 38. Potassium chloride (KCl) (Fisher Scientific, catalog number: P333)
 39. Sodium phosphate monobasic (NaH_2PO_4) (Fisher Scientific, catalog number: S369)
 40. K-gluconate (Sigma-Aldrich, catalog number: P1847)
 41. Phosphocreatine (Sigma-Aldrich, catalog number: P1937)
 42. ATP-Mg (Sigma-Aldrich, catalog number: A9187)
 43. GTP (Sigma-Aldrich, catalog number: G8877)
 44. Poly-L-Lysine (PLL solution) (see Recipes)
 45. Dissection solution (see Recipes)
 46. Papain solution (see Recipes)
 47. DNase/L-cysteine solution (see Recipes)
 48. DNase/Mg solution (see Recipes)
 49. Trypsin inhibitor solution (see Recipes)
 50. Culture medium (NB/B27 medium) (see Recipes)
- B.** For electrophysiology recordings
1. Borosilicate glass capillaries (1.5 OD) (World Precision Instruments, catalog number: 1B150F-4)

2. Artificial cerebrospinal fluid (aCSF, see Recipes)
3. Intracellular solution (see Recipes)

Equipment

- A. For the neuronal culture
 1. Tissue culture hood
 2. P1000 pipet
 3. P200 pipet
 4. Cell culture incubator (Thermo Fisher Scientific, Thermo Scientific™, model: Forma™ Series II 3110, catalog number: 3110)
 5. Vibratome slicer (Leica Biosystems, model: Leica VT1200 S)
 6. Dumont #5 forceps (Fine Science Tools, catalog number: 11251-30)
 7. Dumont #7 forceps (Fine Science Tools, catalog number: 11271-30)
 8. Extra fine bonn scissors (straight) (Fine Science Tools, catalog number: 14084-08)
- B. For the electrophysiology setup
 1. Upright water immersion microscope with fluorescence (Olympus, model: BX51)
 2. Fluorescence filter set for mCherry (TRITC, Chroma Technology, model: 41002c)
 3. CCD camera for fluorescence imaging (Hamamatsu Photonics, model: C8484)
 4. CCD camera for IR imaging (Olympus, model: OLY-150 IR)
 5. Micromanipulator (Luigs & Neumann, model: SM6)
 6. Patch-clamp amplifiers (Dagan, model: BVC-700A)
 7. 18-bit interface card (National Instruments, model: PCI-6289)
 8. Flaming/Brown micropipette puller (Sutter Instrument, model: P-97)
 9. Hemacytometer (Hausser, catalog number: 3110)
- C. For the optogenetic setup
 1. DLP projector (Texas Instruments, model: DLP® LightCrafter™ Evaluation Module)
 2. Aluminum Breadboard, 150 × 300 × 12.7 mm, double density, M6 thread (Thorlabs, model: MB1530/M)

3. Ø1" Achromatic Doublet, SM1-Threaded Mount, $f = 35$ mm, ARC: 400–700 nm (Thorlabs, model: AC254-035-A-ML)
4. Ø2" Achromatic Doublet, SM2-Threaded Mount, $f = 200$ mm, ARC: 400–700 nm (Thorlabs, model: AC508-200-A-ML)
5. Ø2" (Ø50.8 mm) Protected Silver Mirror, 0.47" (12.0 mm) Thick (Thorlabs, model: PF20-03-P01)
6. Kinematic Mount for Ø2" Optics (Thorlabs, model: KM200)
7. U-DP; Dual port intermediate tube (Olympus, model: U-IT140)
8. U-EPA2; Eyepoint adjuster, BX2, Raises eyepoint 30MM (Olympus, model: U-IT101)
9. BX2 filter cube (U-MF2) for Olympus BX2 and IX2 models (Chroma Technology, model: U-MF2, catalog number: 91018)
10. 510 nm beamsplitter (BS, Part Size: $25.5 \times 36 \times 1$ mm) (Chroma Technology, model: T510lpxrxt)
11. 50/50 beamsplitter (BS, Part Size: $25.5 \times 36 \times 1$ mm) (Chroma Technology, model: 21000)
12. EVGA GeForce GT 520 graphic card with a mini HDMI connector (EVGA, model: 01G-P3-1526-KR)
13. Light-to-voltage optical sensors (ams, TAOS, model: TSL13T)

Procedure

A. Primary cultures of neurons expressing channelrhodopsin

Note: Protocol for the culture of primary neurons was originally described by (Brewer et al., 1993) and was further modified by Hilgenberg and Smith (2007) and Barral and Reyes (2016). This procedure requires aseptic conditions for which the Gibco Cell Culture Basics Handbook can provide a good introduction.

1. Coating of the coverslips (prepare the day before culturing under sterile conditions)
 - a. Clean the coverslips in 3 N HCl or concentrated nitric acid (70% wt/wt) for at least 24 h.
 - b. Rinse the coverslips 3 times in distilled H₂O, clean them with 70% ethanol and let dry under the tissue culture hood.
 - c. Coat cleaned coverslips with 0.1 ml of PLL solution (see Recipes) for 12 h.
 - d. Rinse 3 times in distilled H₂O and let dry.

2. Dissociation and culture of primary neurons. Depending on the size of the slicing container, 100 to 150 ml of dissection solution is needed (see Recipes).
 - a. Cut 5–6 600 μm coronal slices of cortex from postnatal (P0–P1) mice in ice-cold dissection solution.
 - b. In the meantime, prepare the papain solution (see Recipes).
 - c. Using forceps, remove the thin membrane covering the brain called the meninges. Dissect out the cortex and cut the cortical layer in small pieces of 0.5–1 mm in length. Transfer the brain pieces in clean dissection solution.
 - d. From this step on, the procedure has to be performed under a tissue culture hood. Filter the papain solution in a Petri dish and incubate the brain pieces in the papain solution for 20 min at 37 °C.
 - e. In the meantime, dilute the DNase/Mg solution (see Recipes) in 10 ml and prepare 4 conical tubes as follow:
 - 1: 4 ml of DNase/Mg solution
 - 2: 2 ml of DNase/Mg solution + 200 μl of trypsin inhibitor solution (see Recipes)
 - 3: 2 ml of DNase/Mg solution + 60 μl of trypsin inhibitor solution
 - 4: 2 ml of DNase/Mg solution + 30 μl of trypsin inhibitor solution
 - f. Using the transfer pipet, transfer the tissue in tube 1 and let the pieces settle down. Try to transfer the brain tissue with the smallest possible volume of solution. Repeat with tubes 2–4.
 - g. Transfer the pieces in about 0.6 ml of equilibrated culture medium (see Recipes) in a low retention tube. Pipet 4–5 times with a P1000 low retention pipet and let the large pieces settle down. Transfer the pellet to a second low retention tube containing about 0.6 ml of culture medium. Pipet 4–5 times with a P200 low retention pipet and pool the two tubes.
 - h. Dilute 10 μl of cell suspension in 10 μl of 0.4% trypan blue and estimate the number of viable cells using a hemacytometer. This is useful to control the final density of neurons.
 - i. Plate $0.2\text{--}3 \times 10^6$ cells on each coverslip, resulting in a density of $\sim 50\text{--}1,000$ cells/ mm^2 at the time of experiment. Let settle for 1 h (in the incubator) for cells to adhere before flooding

the Petri dish with 3 ml of equilibrated culture medium. Maintain cultures at 37 °C and 5% CO₂.

- j.** Exchange one third of the medium every 2–3 days. Cultures can be maintained in these conditions for up to 3 weeks (see Figure 1 for representative results).

3. Expression of channelrhodopsin

- a.** The virus AAV2-hSyn-hChR2(H134R)-mCherry) was produced at 3×10^{12} cfu/ml by the University of North Carolina Vector Core Services using plasmid generously provided by Karl Deisseroth (Stanford University). The virus could also be produced in house using standard protocol (Luo *et al.*, 2007).
- b.** Infect each Petri dish with 1 µl of virus after 3 days *in vitro*.
- c.** Experiments can be performed between 14–21 DIV, when neuronal characteristics and network connectivity were stable and expression of ChR2 was sufficient to enable reliable photostimulation.

B. Building and calibrating the photostimulation setup

The setup for optical stimulation using a video projector was originally described in (Stirman *et al.*, 2012) and further modified in (Barral and Reyes, 2016). The reader will find much details of the procedure to combine a video projector with the microscope in (Stirman *et al.*, 2011; Stirman *et al.*, 2012). The general idea is to focus the image of the projector in the same plane as the sample. Because our microscope was mounted on a movable platform, the photostimulation setup had to be physically connected to the microscope. We describe here the procedure to achieve this but alternative solutions are available and might be more appropriate depending on the setup configuration. Here, we attached a breadboard onto the fluorescent path of the microscope to build the photostimulation light path. We chose the breadboard small enough to limit weight on the microscope but it forced us to use an additional mirror. If possible, try to avoid this mirror by having a longer straight path to project the light coming from the projector.

1. Building the stimulation setup

- a.** Attach the dichroic mirror DM1 to the cube and insert it in the dual port. Place the eyepoint adjuster and the dual port above the fluorescent path and below the camera tube lens of the microscope (Figure 2A).
- b.** Fix the aluminum breadboard on the fluorescent tube and define the optical path. Make use of a reflective mirror (Figure 2A, RM) if necessary.

- c. Place the $\varnothing 2''$ lens that will play the role of the tube lens (Figure 2B, PTL) as close as possible to the dichroic mirror and evaluate the position of the primary image plane.
- d. Place the DLP projector and the $\varnothing 1''$ zoom lens (Figure 2B, ZL) to form an image of the projector exactly at this location. It might be useful to mount the projector tube lens on a movable track in order to finely adjust its position when the microscope objective is changed for example.
- e. It is likely that the collimated beam will be larger than the dual port entry of the microscope. Install the photodiode next to the entry port of the microscope (Figure 3) to collect some of this light from the projector that will be used as a trigger to synchronize photostimulation and electrophysiological recordings. If the light beam is smaller than the entry port, place a 10/90 beamsplitter in-between the projector tube lens (PTL) and the dichroic mirror (DM1) to collect some light from the projector and redirect it to the photodiode system.
- f. For recording of neuronal activity, we used patch-clamp recordings of up to 4 neurons in cell-attached and whole-cell configurations to monitor spikes and post-synaptic potentials, respectively. A useful summary of patch-clamp procedure can be found on [Axolbio](#) webpage. Alternatively, the setup is fully compatible with the use of microelectrode array (MEA) recordings if many neurons need to be recorded. Different commercial MEA systems are available or an open-source DIY system called NeuroRighter can be built for a lower cost.

2. Calibrating the setup

- a. Brightness and pixel size
 - i. Project a full-on pattern and measure the total light flux using a light power meter (Figure 4A).
 - ii. Measure the area of projection using a calibrated camera to deduce the resulting pixel size and the light intensity in mW/mm^2 . With the maximum current to the blue LED, we achieved $15 \text{ mW}/\text{mm}^2$.
- b. Contrast ratio

Here we want to measure the light intensity at the point of stimulation and at other locations to estimate the contrast ratio. This measure is important because it provides an estimate of the light intensity that a non-stimulated neuron receives during stimulation of ChR2 expressing neurons.

- iv. This procedure can be achieved either by placing a mirror at the sample plane. One can also simply look at the reflection of the projector on the glass coverslip of the sample, which is useful if one want to run this calibration during an actual experiment.
 - v. Draw an arbitrary ROI, build the corresponding pattern, project and image it to assess the quality of the calibration.
3. Interfacing the projector and the computer

The projector was controlled by the computer through a USB-based connection. One display mode ('pattern sequence' mode) enabled displaying a limited number of single-bit pattern at a rate up to 4 kHz. With the DLP lightcrafter, only 96 patterns can be stored. In this configuration, the predefined patterns were loaded into the remote memory of the projector and displayed using an external analog trigger. More recent projectors can have additional memory but it will never be sufficient to display patterns at high rate for more than few seconds. Additionally, the time to transfer data to the projector can be prohibitive.

To increase the number of patterns, an alternative method ('HDMI video' mode) was used to stream the images through a high speed HDMI connection. Here, images were sent continuously from the computer to the projector via the graphic card using the HDMI port. We used the 'monochrome 1-bit per pixel' video mode of the projector. This mode allows bypassing any video processing algorithms. A 24-bits RGB image is then considered as a succession of 24 single-bit planes. The resulting time resolution becomes $24 \times 60 \text{ Hz} = 1,440 \text{ Hz}$ for a single-bit image, but without any limitation on the number of displayed patterns.

Synchronizing accurately the display of each frame with the acquisition of electrophysiological signals was a complicated task because it depends on the precise time at which the computer sends an image to the projector. Since the trigger present on the video projector was already used as an input trigger for the 'pattern sequence' mode, we chose to measure directly the light intensity at the output of the projector using a photodiode. Because, the collimated beam was slightly larger than the entry port of the microscope, we could place the photodiode at the edge of the aperture to capture light that did not enter the microscope (see Figure 3). Due to hardware specifications, there is a small 135 μsec delay for every 24-bits plane (*i.e.*, every 16.7 msec) where the image becomes completely dark. We took advantage of this delay, which serves as an accurate trigger for synchronizing photostimulation and electrophysiological recordings offline.

- a. Install the graphic card.
 - b. Install the Graphical User Interface (GUI) provided by Texas Instrument and follow the instructions. Connect the projector using the USB and the HDMI connectors.
 - c. Define the projector as an external display of the computer.
- C. Independent stimulation of neurons using patterned illumination and recording of neuronal activity
1. Transfer the neuronal culture coverslip on the microscope chamber and flood with artificial cerebrospinal fluid (see Recipes).
 2. Take a fluorescent image of the neuronal culture (Figure 5).
 3. Define regions of interest (ROIs) onto neurons that will be activated (Figure 5).
 4. Using a custom-built interface (Labview or Matlab are suitable for this task), build a temporal sequence of images in order to stimulate each ROI with the appropriate temporal pattern (see Table 1 for the key elements of the interface).
 5. Compile a movie from this sequence of images.
Note: Remember to save the movie without any compression since it would destroy the desired pattern sequence.
 6. Pull recording pipets (aim for a pipet resistance of 5–8 M Ω) and fill patch-clamp electrodes with the internal solution (see Recipes).
Note: Here recordings were made in current-clamp configuration, therefore a high pipet resistance was acceptable and had the advantage of being easier to form a seal.
 7. Select a neuron. Patch the neuron and record the electrophysiological signal in cell-attached or whole-cell configuration.
Note: Add QX-314 (Na⁺ channel blocker) to block spikes if excitatory or inhibitory postsynaptic potentials need to be isolated.
 8. Display the stimulation movie through the projector (see Video 1) and record both the electrophysiological signal and the light intensity from the photodiode. Here, signals were filtered at 5 kHz and digitized at 25 kHz using an 18-bit interface card.
 9. Using the data from the photodiode, determine the precise time point where the light stimulation actually started. Realign the electrophysiological to this true starting point and analyze data offline.

Data analysis

- A. Control experiment

To assess the effectiveness of the photostimulation setup in evoking spikes, we stimulated ChR2-positive neurons. We verified that spikes were evoked reliably when a single neuron was stimulated (Figure 6A). Then we stimulated 11 neurons simultaneously and recorded spikes from 2 of them (Figure 6B). Whereas one neuron responded only when stimulated by the light pulses, the other neuron also displayed action potentials in response to the stimulation of neighboring neurons, meaning that it was an integrant component of the recurrent network.

B. Network stimulation

The final experiment consists in stimulating the selected neurons by a controlled stimulus (Figure 7A). Here, the firing rate of the input was fixed at 5 Hz and the spike trains were not correlated. These quantities can be varied at will by the experimenter. Neuronal activity in non-stimulated neurons was monitored either in cell-attached (Figure 7B) or in whole-cell (Figure 7C) configuration to isolate the spiking activity and the postsynaptic inputs, respectively. Cell-attached recordings (Figure 7B) showed that some cells fired robustly (neuron 2) while little activity was evoked in others (neurons 1 and 3). Subsequent whole-cell recordings in the same neurons also showed that large voltage transients evoked reliable action potentials across trials. The amplitudes were much larger than the unitary EPSPs and were mostly likely due to synchronous recurrent activity (Barral and Reyes, 2016).

Notes

1. The protocol to culture primary cortical neurons can work equally on hippocampal neurons.
2. Add 5 μ M Ara-C or 5 μ M 5-Fluoro-2'-deoxyuridine on day 3–5 if glial cells show excessive proliferation. However, it is usually not necessary since Neurobasal culture medium with B-27 supplement does not promote glia cell proliferation and it might increase cell death.
3. We described here a setup using a 460 nm light source to stimulate ChR2-expressing neurons. The setup can easily be modified to use another LED of the projector. Alternatively, the light engine of the projector can be dismantled and an external light source can be used to enlighten the digital micromirror device.

Recipes

1. Poly-L-Lysine (1 \times , 0.1 mg/ml PLL, 20 aliquots of 0.5 ml)
 - a. Make 100 ml of boric acid/sodium tetraborate solution (0.1 M, pH 8.4) Dissolve 0.95 g sodium tetraborate decahydrate in 100 ml ddH₂O Add boric acid to borax solution until desired pH (pH 8.4) is reached (about 0.61 g of boric acid)
 - b. Dissolve 1 mg of PLL in 10 ml boric acid/sodium tetraborate solution

Acknowledgments

Jeremie Barral was supported by a Human Frontier Science Program long-term postdoctoral fellowship (LT000132/2012) and by the Bettencourt Schueller Foundation. Alex Reyes was supported by grants from the National Institutes of Health (DC005787-01A1). This protocol was adapted from procedures published in Barral and Reyes (2016).

References

- Anselmi, F., Banerjee, A., Albeanu, DF. Patterned photostimulation in the brain. In: Douglass, AD., editor. *New Techniques in Systems Neuroscience* Springer. 2015. p. 235-270.
- Barral J, Reyes A. Synaptic scaling rule preserves excitatory-inhibitory balance and salient neuronal network dynamics. *Nat Neurosci.* 2016; 19(12):1690–1696. [PubMed: 27749827]
- Boyden ES, Zhang F, Bamberg E, Nagel G, Deisseroth K. Millisecond-timescale, genetically targeted optical control of neural activity. *Nat Neurosci.* 2005; 8(9):1263–1268. [PubMed: 16116447]
- Brewer GJ, Torricelli JR, Evege EK, Price PJ. Optimized survival of hippocampal neurons in B27-supplemented neurobasal, a new serum-free medium combination. *J Neurosci Res.* 1993; 35(5): 567–576. [PubMed: 8377226]
- Farah N, Reutsky I, Shoham S. Patterned optical activation of retinal ganglion cells. *Conf Proc IEEE Eng Med Biol Soc.* 2007; 2007:6368–6370. [PubMed: 18003478]
- Guo ZV, Hart AC, Ramanathan S. Optical interrogation of neural circuits in *Caenorhabditis elegans*. *Nat Methods.* 2009; 6(12):891–896. [PubMed: 19898486]
- Hilgenberg LG, Smith MA. Preparation of dissociated mouse cortical neuron cultures. *J Vis Exp.* 2007; (10):562. [PubMed: 18989405]
- Luo J, Deng ZL, Luo X, Tang N, Song WX, Chen J, Sharff KA, Luu HH, Haydon RC, Kinzler KW, Vogelstein B, He TC. A protocol for rapid generation of recombinant adenoviruses using the AdEasy system. *Nat Protoc.* 2007; 2(5):1236–1247. [PubMed: 17546019]
- Nawrot MP, Schnepel P, Aertsen A, Boucsein C. Precisely timed signal transmission in neocortical networks with reliable intermediate-range projections. *Front Neural Circuits.* 2009; 3:1. [PubMed: 19225575]
- Papagiakoumou E, Anselmi F, Begue A, de Sars V, Gluckstad J, Isacoff EY, Emiliani V. Scanless two-photon excitation of channelrhodopsin-2. *Nat Methods.* 2010; 7(10):848–854. [PubMed: 20852649]
- Shoham S, O'Connor DH, Sarkisov DV, Wang SS. Rapid neurotransmitter uncaging in spatially defined patterns. *Nat Methods.* 2005; 2(11):837–843. [PubMed: 16278654]
- Stirman JN, Crane MM, Husson SJ, Gottschalk A, Lu H. A multispectral optical illumination system with precise spatiotemporal control for the manipulation of optogenetic reagents. *Nat Protoc.* 2012; 7(2):207–220. [PubMed: 22240583]
- Stirman JN, Crane MM, Husson SJ, Wabnig S, Schultheis C, Gottschalk A, Lu H. Real-time multimodal optical control of neurons and muscles in freely behaving *Caenorhabditis elegans*. *Nat Methods.* 2011; 8(2):153–158. [PubMed: 21240278]




Figure 1. Culture of primary cortical neurons at 14 days *in vitro*

Representative examples of a healthy culture where cell bodies are well defined (A) and a culture where the coverslip coating was not achieved properly which resulted in cell aggregates and debris (B).




Figure 2. Photostimulation setup

A. Picture of the photostimulation setup; B. Schematic of the light path. DLP, Digital Light Processing projector; ZL, zoom lens; PTL, projector tube lens; EPA, eyepoint adjuster; CTL, camera tube lens; DM, dichroic mirror; RM, reflective mirror.




Figure 3. Photodiodes system

The photodiodes are placed next to the entry port of the microscope and collect a fraction of the beam coming from the projector (dashed line).




Figure 4. Calibration of the photostimulation setup

Images sent to the projector (left) and resulting images from the camera (right) when a full-on (A), ANSI checkerboard (B), or a 4×4 white dots (C) pattern were projected.




Figure 5. Pattern stimulation of the neuronal culture

36 regions of interest (ROIs) are drawn around stimulated neurons (blue squares). 3 patch-clamp electrodes record spikes and/or subthreshold membrane potential of non-stimulated neurons.




Figure 6. Pattern stimulation efficacy in evoking spikes in ChR2-expressing neurons

In the following schematics, neurons expressing ChR2 are shown in light red. Recorded neurons are designated by the recording pipets. Light-stimulated neurons have colored contours. A. A train of action potentials (black dots) was applied to two different neurons expressing ChR2 (blue and orange) independently (top and bottom, respectively). Several trials were realized to confirm the faithful stimulation of neurons of interest. Colored dots denote recorded action potential and grey lines the light flashes. No action potentials were observed in the non-stimulated neuron. B. Then, we stimulated 11 neurons in the network. We generated 11 spike trains (top) that were applied to neurons expressing ChR2, including the two previous neurons. Here, we used correlated spike trains for the stimulus (the average correlation between spike trains was $C = 0.5$). The blue and the orange neurons were recorded simultaneously.




Figure 7. Independent stimulation of the neuronal network

A. Spike trains which are used for stimulating the 36 selected neurons in Figure 5. Each line is a train of light pulses applied to a given ROI. B. Raster plot of three neurons recorded simultaneously. Each line represents a given repetition of the same stimulus (15 repetitions). C. Corresponding membrane potential of Cell 2. Each repetition (7 repetitions) is drawn as a grey line and the average is displayed in red. The patch pipet contained the Na⁺ channel blocker QX-314 to block spikes and isolate postsynaptic potentials.

**Video 1. Pattern stimulation of the neuronal culture**

36 regions of interest (ROIs) are drawn around stimulated neurons (blue squares). Each stimulated neuron is activated by a train of light pulses (the actual stimulus is shown in blue on the top left corner). Three neurons are recorded simultaneously in cell-attached mode (yellow raster plots of evoked spikes in each corner, 15 repetitions of the stimulus) and then in whole-cell configuration (yellow trace of membrane potential, average of 7 repetitions).

Table 1

Key elements that the custom-built interface must be able to achieve

	Goal	Procedure
1	Calibrate the DLP and the camera	Match camera pixels to DLP pixels
2	Determine ChR2-positive neurons	Record fluorescent image with the camera
3	Select neurons to be stimulated	Define ROIs onto neurons to stimulate
4	Define the sequence of activation	Generate Poisson-like spike trains
	Build the corresponding patterns	Assign each spike train to each neuron
5	Save the sequence as a movie	Combine all patterns
6	Stimulate neurons	Display the movie
7	Acquire data	Record electrophysiological and photodiode data

Author Manuscript

Author Manuscript

Author Manuscript

Author Manuscript