Sample	Detector	СТ	Mean	DCT	DDCT	Ratio sample/ref
Hyperdiploid	ZCWCC2	28.22	28.29	13.915	-2.72	6.6
	18S	14.42	14.375			
Hyperdiploid	ZCWCC2	26.74	26.775	11.985	-4.65	25.1
	18S	14.74	14.79			
Other	ZCWCC2	35.47	35.44667	17.95167	1.316667	0.4
	18S	17.62	17.495			
Other, reference	ZCWCC2	31.62	31.64	16.635	0	1
	18 <b>S</b>	15.01	15.005			

 Table 7. Results from the Taq-Man assay for the ZWCC2 gene for the leukemias; two samples were positive for the hyperdiploid, and two negative