

Table 1. QRT-PCR Conditions

Target	5' primer	3' Primer	Annealing temp, °C
GLUT4	5'-cagctctcaggcatcaat-3'	5'-tctactaagagcaccgag-3'	53
CHOP10	5'-aatctggagagcgagggtt-3'	5'-tcactactcttgaccctgcg-3'	57
C/EBP α	5'-tggacaagaacagcaacgag-3'	5'-tcaactggtaactccagcac-3'	57
BiP	5'-cagatcttctccacggcttc-3'	5'-ttcagctgtcactcggagaa-3'	55
18S	5'-cgaccaaaggaaccataact-3'	5'-tggttgatcctgccagtag-3'	57