



Fig. S1. Aggregation of β -amyloid peptide (A β 1-42) in the presence of *cyt c*/TOCL/ H_2O_2 analyzed by SDS-PAGE with subsequent Western blotting (staining with anti- β -amyloid peptide antibody). Accumulation of covalently cross-linked hetero-oligomers was observed when A β and *cyt c* plus TOCL were incubated in the presence of H_2O_2 (lane 4, arrow). Conditions: A β (2 μ M) was incubated with *cyt c* (2 μ M) and TOCL/DOPC liposomes (50 μ M, TOCL/*cyt c* ratio 25:1) for 60 min at 37°C. H_2O_2 (50 μ M) was added every 15 min. Excess H_2O_2 was removed by the addition of catalase (1400 units/ml catalase). A typical gel representative of three experiments is shown.