



Supplemental Fig. 2. CXCR3 expression on endothelial tubes

To verify that CXCR3 is expressed on newly formed tubes, HMEC-1 cells were plated on 8 well chamber cover slips (Nalgene) coated with growth factor reduced Matrigel (35,000 cells/well). Cells were incubated for 48 hrs in 2.0% dialyzed FBS-MCDB131 media. The cells were fixed with 2% paraformaldehyde-PBS then immuno stained with mouse anti-CXCR3^{FITC} (Green) antibody (R&D Systems, 1:300 dilution). Cells were also incubated with Draq5 (nuclei, Blue). The tubes were imaged on an Olympus Fluoview confocal microscope using a LucPlanFLN 20x/0.45 objective.