



Rxn No.	NEL2-link 5' (928 bp)	NEL2-1 (2392 bp)	NEL2-2 (3037 bp)	NEL2-3 (4089 bp)	NEL2-link 3' (1914 bp)
JGS1651	+	+	+	+	+
JGS5252	+	+	+	+	+
JGS5621	+	+	+	+	+
JGS1521	+	+	+	+	+
JGS4143	+	+	+	+	+
JGS4140	+	+	+	+	+
CP4	+	+	+	+	+
CP1	+	+	+	+	+
CP2	+	+	+	+	+
CP3	+	+	+	+	+
CP6	+	+	+	+	+
JGS1473	-	-	-	-	-