



Rxn No.	NEL3-1 (2.9 kb)	NEL3-2 (940 bp)	NEL3-3 (2.9 kb)	NEL3-4 (1.9 kb)
JGS1651	+	+	+	+
JGS5252	+	+	+	+
JGS5621	+	+	+	+
JGS1521	+	+	+	+
JGS4143	+	+	+	+
JGS4140	+	+	+	+
CP4	+	+	+	+
CP1	+	+	+	+
CP2	+	+	+	+
CP3	+	+	+	+
CP6	+	+	+	+
JGS1473	-	-	-	-