

	50% titer				endpoint titer (95%)			
	VarO	R29	PF13	PFL	VarO	R29	PF13	PFL
<b>Anti-VarO</b>	<b>409 600</b>	4 000	10 00	1 000	<b>6 553 600</b>	204 800	409 600	6 400
<b>Anti-R29</b>	4 000	<b>204 800</b>	5 000	2 000	102 400	<b>3 276 800</b>	40 000	12 800
<b>Anti-PF13</b>	1 000	2 000	<b>204 800</b>	2 500	204 000	75 000	<b>3 276 800</b>	25 600

	50% titer				endpoint titer (95%)			
	VarO	R29	PF13	PFL	VarO	R29	PF13	PFL
<b>Anti-VarO</b>	<b>409 600</b>	100-fold	40-fold	400-fold	<b>6 553 600</b>	30-fold	16-fold	1000-fold
<b>Anti-R29</b>	50-fold	<b>204 800</b>	40-fold	100-fold	30-fold	<b>3 276 800</b>	80-fold	250-fold
<b>Anti-PF13</b>	20-fold	100-fold	<b>204 800</b>	80-fold	16-fold	45-fold	<b>3 276 800</b>	128-fold