Table S4. Proportions of nascent proinsulin monomers and non-monomers precipitated by insulin (Ins) or C-peptide (Cp) antisera from mouse islets labeled for 30 minutes

Percentage	Proinsulin State	Ins IP	Cp IP
Mean	Monomers	58.2	39.7
Mean	Non-monomers	41.8	60.3
SD	Monomers	4.2	7.6
SD	Non-monomers	4.2	7.6
P (Monomers; Ins IP vs. Cp IP)			<0.005
P (Non-monomers; Ins IP vs. Cp IP)			<0.005

(Shown in Figure 2B)