

### **RT-PCR primers**

p63Alpha/Beta: For	CAGCACCAGCACTTACTTCAGA
p63 alpha Rev: Rev	CCACCTCGGTAAGAACTGA
p63 Beta: Rev	TCAGACTTGCCAGATCCTGA
p63 Gamma: For	ACAGGAAGGCGGATGAAGATAGCA
p63 Gamma: Rev	GGCTCCACAAGCTCATTCT
p130: For	CAGGCCCGGAGCCAGGTGTA
p130: Rev	GGAGCACTGCTGGGCTGTGG
Rb: For	CTCCACAGCTCGCTGGCTCC
Rb: Rev	CTCAAGCCTGACGAGAGGCAGGT
p107: For	CCCCGCACAAGAATGGGTCAGG
p107: Rev	ACAGACGCGTTTGGCAGGGG
Skp2: For	TGCGGACGCTATGCACAGGAAG
Skp2: Rev	GCGTGGCCAGGTTTGA
RPLPO: For	ATCAACGGGTACAAACGAGTC
RPLPO: Rev	CAGATGGATCAGCCAAGAAGG

### **ChIP primers**

Skp2 -234: FOR-	ATCCGTTCTACAGTCCAGC
Skp2 -234: REV	TACCCGCCTTCGAGATAC
Skp2 -1254: FOR	TTCCTCAAGAGGCAGTTA
Skp2 -1254: REV	CGCATATTAGAATATTCAGTG
Skp2 intron-1: For	GGAAGCCTGGGCCTCCAGTG 3'
Skp2 intron-1: Rev	GGAAGCCTGGGCCTCCAGTG 3'