Elective Selection of Proline-requiring Mutants

YUKINORI HIROTA, MANABU INUZUKA, AND MUNEMITSU TOMOEDA

Department of Biology, School of Science, Osaka University, Osaka, and Faculty of Pharmaceutical Sciences, University of Kanazawa, Kanazawa, Japan

Received for publication 3 January 1966

We recently discovered a method whereby proline-requiring (pro) mutants of some strains of Escherichia coli K-12 can be selected at high frequency (Table 1). The method seems to involve the general principle of differential growth rates. Only pro E. coli K-12 mutant strains grew in Penassay Broth (Difco) containing 1 to 3 μ g/ml of 4-nitropyridine-N-oxide (4NPO), 10 μ g/ml of 4-nitrosopyridine-N-oxide, or 50 μ g/ml of 4-hydroxyl-aminopyridine-N-oxide, whereas strains not requiring proline (pro^+) were inhibited.

The method can be illustrated by the results obtained with a pro+ E. coli K-12 strain, W3630. A single colony of this strain was inoculated into 10 ml of Penassay Broth and incubated at 37 C overnight, during which time growth usually reached a level of about 5×10^8 cells per milliliter. The cell suspension was diluted 10⁻³, and 0.2 ml was added to 100 ml of Penassay Broth along with 0.2 to 0.5 ml of a solution of 1 mg of 4NPO/ml. This mixture was incubated on a reciprocal shaker at 37 C for 20 to 26 hr, and 0.1 ml of the culture was diluted and plated on eosin-methylene blueglucose-agar medium (J. Lederberg, Methods Med. Res. 3:55, 1950). Colonies appearing on this medium then were replica-plated onto Davis's minimal (DM) agar with or without 40 μ g/ml of L-proline. The percentage of pro mutants present in 4NPO-treated cultures is shown in Table 1. Pretreatment of bacteria with a mutagen, such as ultraviolet light or N-methyl-N'-nitro-N-nitrosoguanidine, increased the frequency of appearance of pro mutants.

A preliminary analysis of the method has been made with the following results: pro mutants are the only ones selected by this method at a detectable frequency; pro mutants of independent origin grow on DM agar containing the metabolic precursors of proline, glutamic-γ-semialdehyde or

 $\Delta 1$ -pyrroline-5-carboxylic acid, but they do not grow on DM agar containing glutamic acid (presumably the mutants have a metabolic block between glutamic acid and glutamic- γ -semialdehyde); the mutants produced are stable; conjugation experiments indicate that the *pro* marker is located on the chromosome between the arabinose (*ara*) and lactose (*lac*) loci, close to *lac*. Re-

Table 1. Selection of proline-requiring mutants from Escherichia coli K-12 W3630

Concn of 4NPO	No. of pro mutants*				
0	0	0	0	0	0
1.0	18.5	8.9	38.5	0.9	17.9
2.0	59.3	53.3	100		
3.0	93.9	1.7	0	0	

* Each culture contained about 10^4 to 10^8 cells per milliliter; the number of colonies tested in each experiment was about 200. Results are expressed as: (number of pro mutants)/(number of colonies tested) \times 100.

construction experiments with a mixture of pro lac⁺ and pro⁺ lac⁻ cells showed that the appearance of pro mutants was not due to mutation induced by 4NPO, but rather to selective growth in the media.

We thank R. Sato for fruitful discussions, and Y. Izumi for a gift of glutamic- γ -semialdehyde dimethyl acetal. We also thank S. Kanahara, A. Oiwa, and N. Hiyoshi of the Faculty of Pharmaceutical Sciences, Kanazawa University, for preparation of 4NPO and its derivatives.

This investigation was supported by Public Health Service research grant GM 08293 from the Division of General Medical Sciences.