

VISTA Enhancer Browser

enhancer ID	fragment length (bp)	Gene	fwd primer (5' - 3')	rev primer (5'-3')	Restriction enzymes used to release the insert for pronuclear injection
hs611	854	Cux2	611fwd: AAAC <del>T</del> CGAGCTCTTCTCTCTGTTCCTGTGGGTC	611rev: AAAGTCGACCCGACCGTTCCCTCTTCCCT	XhoI / NotI
hs112	1873	Dmrt3	112.fw3: GACACTGATAGGGTCAGTGC	112.rev3: AACTAGTTTCTCTGGGAAAGAGAG	NotI
			112fw2: GGAAACGGAGTTCATGTGGT	112rev: AAAGTCGACTACCGCCAGCCAAGGAAGAGAGCA	NotI
hs111	438	Gli3	111fwd: AAAGGTACCGAGCAATTGCGGAGTGCAGGCC	111rev: AAAAAGCTTCATTACCTCTCTCCAACCAGC	KpnI / NotI
hs1450	823	Nfia	NFIA_fw: GCCTGTTTCCTGCCGACCTGA	NFIA_rev: CTCCTTGGA <del>A</del> CTCACAGTCAGGGA	
hs1006	1068	Wnt8b	Wnt8b_F1: GACCCCTTCCTCCTTGT	Wnt8b_R1: GGTTTCCCACACCACTTC	