

Supplementary Table 4. Primers for cloning the ABCA1 promoter

Target vector	Forward primer (5' → 3')	Reverse primer (5' → 3')	Amplicon size (bp)
<i>pCR2.1</i>	TGT TCC TCT CGG GTC CTC T	ACT CAC TCT CGC TCG CAA TTA	715
<i>pGL4.10</i>	GCT TTG TTC CTC <u>TGA</u> <u>TAT</u> CCT CTG A	GCA GAA TTA AGC <u>TTA</u> CTC ACT CTC G	733

The engineered EcoRV and HindIII restriction sites in the forward and reverse pGL4.10 cloning primers respectively are underlined.