

Table S2. Primers for bisulfite sequencing (bisulfite-PCR)

	Forward primer	Reverse primer
<i>Atp8b1</i> top	5'-AAGTTAGTTGTTGTATTAGAAGATTATTAT-3'	5'-CTAACCTTATCTTTACATTAACAAAAA-3'
<i>Atp8b1</i> bottom	5'-TTGGTTTTGTTTTGTATTGAAGAG-3'	5'-ACATCAAAAAATTATCACAAAAATC-3'
<i>Igf1r</i> top	5'-TGATTGGATATTTGATAAGATTATTTTGA-3'	5'-AAAAAATACAACATCCAACCAAAAA-3'
<i>Igf1r</i> bottom	5'-ATTTTGGTGATTGAAAAGTTTATTG-3'	5'-CCAATTCATCATACTAATACTATATCTT-3'
<i>Igf2r</i>	5'-GGGAATTGAGGTAAGTTAGGGTTTT-3'	5'-TCTTATAACCCAAAAATCTTCACCCTAAC-3'
<i>H19</i>	5'-GTTTGAGGAGTTTTAAGGTAGAAGGGGATT-3'	5'-TAACAAACTCATAAATCACTCAAACATAACATTC-3'