
Supplemental Figure S2



BAFF-induced processing of NF-κB2 p100 to p52 in Bcl-3-/- B cells.

B cells isolated from WT and KO mice were stimulated with BAFF (500 ng/ml) (Peprotech) for the times indicated. Western blotting was performed to examine processing of p100 to p52. Freshly isolated cells contain processed p52 due to exposure to BAFF *in vivo*; most of the *in vivo* processed p52 disappears after an overnight incubation and will reappear only upon addition of BAFF to the culture. A second independent experiment yielded similar results.