

Table S1 : MS Identification of protein spots of 2D-PAGE

protein ID	mass weight	pI	peptide count	Mascot protein score	Mascot protein score C.I. %
CysC	22503	5,25	4	93,7	100
GamA	27271	5,81	6	80,8	100
LiaH	25682	6,2	10	261	100
LiaR	23110	5,13	12	158	100
NadE	30376	5,07	16	456	100
NrfA	28302	5,78	12	224	100
PbpC	74360	6,24	19	373	100
PspA	25125	5,87	12	85,1	100
RpsB	27949	6,27	11	207	100
YceC	21809	5,46	11	110	100
YceH	41646	5,9	20	303	100
YpuA	31275	4,51	38	20679*	-
YqiG	40779	5,34	13	277	100
YtrE	25444	5,99	13	198	100
YvcR	29020	6,32	10	261	100
YvIB	41055	5,5	15	123	100

*identified with nanoUPLC-ESI-MS/MS with a Synapt G2S mass spectrometer. The PLGS score is given.