

**SUPPLEMENTARY FIG. S2.** RNase refolding over time with hBCAT +/- protein disulphide isomerase (PDI). PDI and each respective hBCAT isoform were incubated with rdRNase ( $30 \,\mu\text{M}$ ), and aliquots were removed at 24 h to determine the effect on RNase folding. Change in absorbance at 295 nm over time for PDI ( $\blacksquare$ ), hBCATc ( $\bigcirc$ ), hBCATm ( $\blacktriangledown$ ), PDI+hBCATc ( $\triangle$ ), and PDI+hBCATm ( $\blacksquare$ ). Graph is representative of three independent runs.