

Figure S1




Fig S1. Effects of heteronemin on cell cycle distribution in A498 cells. Cells were incubated with vehicle or various concentrations of heteronemin for 24 h or indicated time and detected cell cycle distribution by flow cytometry.

Figure S2




Fig S2. Effects of SB203580 on heteronemin-induced LC3 conversion in A498 cells. Cells were incubated with vehicle, heteronemin 3 μ M, SB203580 25 μ M and combination of heteronemin and SB203580 for 24 h and detected protein expression by western blotting.