

**S3 Table. Details of closed-loop prosthetic hand control.**

Subject	Decoding of movement type		Detection of movement onset				<i>p</i> -value
	# of correct movements	# of total movements	# of TP	# of FP	# of FN	# of TN	
1	10	12	15	6	2	11	0.002
2	13	18	17	12	3	8	0.078
3	10	13	14	5	7	12	0.024
4	11	21	16	10	9	10	0.261
5	6	10	10	6	9	9	0.350
6	11	14	14	6	4	12	0.009

TP: true positive, TN: true negative, FP: false positive, FN: false negative for onset detection; *p*-value: *p*-value of one-sided Fisher's exact test.