

S6 Table. SmartCube results for the Cntnap2 knockout model.

Cntnap2					
SmartCube	Measure		Genotype	Mean	SE
	Repetitive Behavior	Grooming Time	WT	193.9	21.6
			KO	185.6	20.1
		Digging Time (s)	WT	0.3	0.04
			KO	0.2	0.03
	Exploration	Latency to Approach Aversive Stimulus	WT	36.0	7.4
			KO	13.5	2.8
		Sniffing (frequency)	WT	446.9	29.3
			KO	396.1	25.4
		Unsupported Rearing Time (s)	WT	72.5	7.1
			KO	65.2	9.1
	Activity	Freezing Time (s)	WT	1.8	0.3
			KO	0.4	0.1
		Abrupt Movement Frequency	WT	6.8	1.3
			KO	30.7	5.4
		Velocity	WT	3.3	0.1
			KO	3.8	0.1
		Locomotion Bursts (frequency)	WT	236.6	14.9
			KO	358.5	21.7
	Short Steps (frequency)	WT	42.4	3.4	
KO		62.8	2.9		