

**Table S1. Primers used for qRT-PCR**

Primers used for qRT-PCR	5' -> 3'
Rp49 Forward	GACGCTTCAAGGGACAGTATCTG
Rp49 Reverse	AAACGCGGTTCTGCATGAG
Drosomylin (Drs) Fw	CGTGAGAACCTTTTCCAATATGATG
Drosomylin (Drs) Rv	TCCCAGGACCACCAGCAT
Drosocin (DroA) Fw	CACCCATGGCAAAAACGC
Drosocin (DroA) Rv	TGAAGTTCACCATCGTTTTCTG
Diptericin (Dipt) Fw	GCTGCGCAATCGCTTCTACT
Diptericin (Dipt) Rv	TGGTGGAGTGGGCTTCATG
Attacin (AttC) Fw	TGGGCTACAACAATCATGGA
Attacin (AttC) Rv	GCGTATGGGTTTTGGTCAGT
Cecropin (Cec) Fw	ACGCGTTGGTCAGCACACT
Cecropin (Cec) Rv	ACATTGGCGGCTTGTTGAG
Metchnikowin (Mtk) Fw	CGTCACCAGGGACCCATTT
Metchnikowin (Mtk) Rv	CCGGTCTTGTTGGTTAGGA