



Supplementary Fig. 4. HIVxCD3 DARTs reduce cell-associated HIV p24 protein, HIV RNA, and integrated HIV DNA. PGT121xCD3 + 7B2xCD3 HIV DARTs (200 pM each) were added to unstimulated primary CD4 T cells 24 hours post spinfoction with HIV BaL. After 5 days of culture, cells were analyzed: by FACS for intracellular p24 protein expression, by qRT-PCR (COBAS) for vRNA, and by qPCR for total vDNA. The average of data from 2 donors are depicted.