

Supplementary Table S1 Semen and prepared sperm analysis for Patient 1, Patient 2 and Donor 1.

	Semen		Prepared sample		
	Count (million/ml)	Motile (%)	Count (million/ml)	Motile (%)	Hyperactivation (%)
Patient 1	66 (54.65)	60 (16.5)	41 (33.85)	92 (5)	11 (2)
Patient 2	80 (2)	55 (2.5)	20 (1.5)	87 (4.5)	28 (1.25)
Donor 1	109 (8.85)	69 (16)	36 (2.4)	91 (1)	22 (1)

Data are mean (\pm SD).