

Alias	Cell Lines	Reprogramming method	Reprogrammed cells	Provided by	DMD gene mutations
hESCs 1	SA001	-	-	Cellartis AB	-
hESCs 2	H9	-	-	Wicell research institute	-
hiPSCs 1	VAX1024	Amphotropic retroviruses	Myoblasts	Celogos	-
hiPSCs 2	M00180	Ecotropic retroviruses	Myoblasts	Cochin Hospital-Cochin Institute	-
hiPSCs 3	M00194	Ecotropic retroviruses	Myoblasts	Cochin Hospital-Cochin Institute	-
hiPSCs 4	M00398	Ecotropic retroviruses	Myoblasts	Cochin Hospital-Cochin Institute	-
DMD hiPSCs 1	M00197	Ecotropic retroviruses	Myoblasts	Cochin Hospital-Cochin Institute	In-frame Duplication exons 3-26
DMD hiPSCs 2	M00202	Ecotropic retroviruses	Myoblasts	Cochin Hospital-Cochin Institute	Out-of-frame Deletion exons 8-43
DMD hiPSCs 3	M00418	Ecotropic retroviruses	Myoblasts	Cochin Hospital-Cochin Institute	Stop exon 7 c