## **Guide to Supplementary Information**

- 1. Supplementary Tables 1-6
- 2. Supplementary Videos 1-2 and legends
  - 1. Video 1 (QuickTime 17.1M)
  - 2. Video 2 (QuickTime 14.3M)

These bright-field movies show dislocation of the lens placode and formation of ectopic lens in the *hir* mutant embryo (Video 2) as compared with WT embryo (Video 1).

- 3. Supplementary Discussion 1: ARHGAP18 related genes
- 4. Supplementary Discussion 2: Why does the *hir* medaka mutant exhibits a unique flattened phenotype rather than a cell proliferation defect?
- 5. Supplementary Figure 1: Full-scan of Western blots