



Supplemental Figure S8. Visualization of certificates in IGV. (A) Isoform Iso_743 has been defined and tracked by SpliceHunter based on reads shown for different time points. Iso_743 has been assigned to gene SPAC1296.03c(sxa2) (negative strand) based on overlap. SpliceHunter reports a novel acceptor (Definition 6) and a novel exon (Definition 8) for this isoform. (B) SpliceHunter clustered reads shown for different time points by the common intron to define novel TU NG_I_35 which does not overlap any annotated gene.