

Supplementary Table S1. The used bacteria and plasmids

Bacterium and plasmid	Genotype, property or sequence	Source or reference
Bacterial strains		
<i>E. coli</i> BL21 (DE3)	F ⁻ ompT hsdS _B (r _B ⁻ m _B ⁻) gal (λ1857 ind1 Sam7 nin5 lacUV5 T7gene1) dcm	Novagen
Ec (Cppdo2-Cpsqr)	<i>E. coli</i> BL21(DE3)/pBBR1MCS-2::Cppdo2-Cpsqr	(Xin et al 2016)
Ec (Cppdo2)	<i>E. coli</i> BL21(DE3)/pBBR1MCS-2::Cppdo2	This study
Ec (Cpsqr)	<i>E. coli</i> BL21(DE3)/pBBR1MCS-2::Cpsqr	This study
Ec(Papdo-Pasqr ₂)	<i>E. coli</i> BL21(DE3)/pBBR1MCS-5::Papdo-Pasqr ₂	This study
Ec(Zppdo-Pasqr ₂)	<i>E. coli</i> BL21(DE3)/pBBR1MCS-5::Zppdo-Pasqr ₂	This study
Ec(Sapdo-Pasqr ₂)	<i>E. coli</i> BL21(DE3)/pBBR1MCS-5::Sapdo-Pasqr ₂	This study
Ec(Bcpdo-Pasqr ₂)	<i>E. coli</i> BL21(DE3)/pBBR1MCS-5::Bcpdo-Pasqr ₂	This study
Ec(Zppdo)	<i>E. coli</i> BL21(DE3)/pBBR1MCS-5::Zppdo	This study
Ec(Sapdo)	<i>E. coli</i> BL21(DE3)/pBBR1MCS-5::Sapdo	This study
Ec(Bcpdo)	<i>E. coli</i> BL21(DE3)/pBBR1MCS-5::Bcpdo	This study
Ec(Pasqr ₂)	<i>E. coli</i> BL21(DE3)/pBBR1MCS-5::Pasqr ₂	This study
Ec(pMCS5)	<i>E. coli</i> BL21(DE3)/pBBR1MCS-5	This study
<i>E. coli</i> HB101/pRK2013	For vector mobilization	(Zhang et al 2013)
<i>Corynebacterium vitaeruminis</i> DSM20294	Wild type, Cv as abbreviation	DSMZ
<i>Cupriavidus pinatubonensis</i> JMP134	Wild type, Cp as abbreviation	ATCC 23344
<i>C. pinatubonensis</i> 2K	The pdo2 and sqr mutant of <i>C. pinatubonensis</i> JMP134	This study
<i>Pseudomonas aeruginosa</i> PAO1	Wild type, Pa as abbreviation	ATCC 15692
<i>P. aeruginosa</i> 3K	The pdo1, pdo2 and sqr mutant of <i>P. aeruginosa</i> PAO1	This study
<i>Agrobacterium tumefaciens</i> C58	Wild type, At as abbreviation	ATCC 33970
<i>Bacillus cereus</i> ATCC 10876	Wild type	ATCC 10876
<i>Bacillus subtilis</i> 168	Wild type, Bs as abbreviation	ATCC 23857

<i>Burkholderia cepacia</i> ATCC 25416	Wild type	ATCC 25416
<i>Corynebacterium glutamicum</i> RES167	Wild type	(Tauch et al 2002)
<i>Enterobacter cloacae</i>	Wild type	Yuzhen Zhang
<i>Enterococcus faecalis</i> CMCC32223	Wild type	ATCC 29212
<i>Gluconobacter oxydans</i> 621H	Wild type	DSM 2343
<i>Klebsiella pneumonia</i>	Wild type, Kp as abbreviation	Yuzhen Zhang
<i>Pseudomonas putida</i> S16	Wild type	ATCC BAA2546
<i>Serratiamar cescens</i> ATCC13880	Wild type	ATCC 13880
<i>Sinorhizobium meliloti</i> 1021	Wild type	(Galibert et al 2001)
<i>Sphingobium chlorophenolicum</i> L-1	Wild type	ATCC 53874
<i>Serratia fonticola</i> DSM 4576	Wild type, Sf as abbreviation	DSMZ
<i>Staphylococcus aureus</i> ATCC 6538P	Wild type	ATCC 6538P
<i>Staphylococcus sciuri</i>	Wild type	Yuzhen Zhang
<i>Streptomyces coelicolor</i> M145	Wild type	ATCC BAA471
<i>Streptomyces griseus</i>	Wild type	Yuzhen Zhang
<i>Zunongwangia profunda</i> SM-A87	Wild type	DSM 18752
<i>Zymomonas mobilis</i> ATCC 31821	Wild type	ATCC 31821
Plasmids		
pBBR1MCS-5	Gm ^r , Plac	(Kovach et al 1995)
pBBR1MCS-2::Cppdo2-Cpsqr	Cppdo2 and Cpsqr in pBBR1MCS-2	(Xin et al 2016)
pBBR1MCS-2::Cppdo2	Cppdo2 in pBBR1MCS-2, modified from pBBR1MCS-2::Cppdo2-Cpsqr	This study
pBBR1MCS-2::Cpsqr	Cpsqr in pBBR1MCS-2, modified from pBBR1MCS-2::Cppdo2-Cpsqr	This study
pBBR1MCS-5-Papdo-Pasqr2.	Papdo and Pasqr2 in pBBR1MCS-5, control by IPTG induced lac promoter	This study
pBBR1MCS-5-Zppdo-Pasqr	Zppdo and Pasqr2 in pBBR1MCS-5, control by IPTG induced lac promoter	This study
pBBR1MCS-5-Sapdo-Pasqr	Sapdo and Pasqr2 in pBBR1MCS-5, control by IPTG induced lac promoter	This study
pBBR1MCS-5-Bcpdo-Pasqr	Bcpdo and Pasqr2 in pBBR1MCS-5, control by IPTG induced lac promoter	This study

pBBR1MCS-5:: <i>Zppdo</i>	<i>Zppdo</i> in pBBR1MCS-5, modified from pBBR1MCS-5:: <i>Zppdo-Pasqr2</i>	This study
pBBR1MCS-5:: <i>Sapdo</i>	<i>Sapdo</i> in pBBR1MCS-5, modified from pBBR1MCS-5:: <i>Sapdo-Pasqr2</i>	This study
pBBR1MCS-5:: <i>Bcpdo</i>	<i>Bcpdo</i> in pBBR1MCS-5, modified from pBBR1MCS-5:: <i>Bcpdo-Pasqr2</i>	This study
pBBR1MCS-5:: <i>Pasqr2</i>	<i>Pasqr2</i> in pBBR1MCS-5, modified from pBBR1MCS-5:: <i>Zppdo-Pasqr2</i>	This study

Reference

- Galibert F, Finan TM, Long SR, Puhler A, Abola P, Ampe F *et al* (2001). The composite genome of the legume symbiont *Sinorhizobium meliloti*. *Science* **293**: 668-672.
- Kovach ME, Elzer PH, Hill DS, Robertson GT, Farris MA, Roop RM, 2nd *et al* (1995). Four new derivatives of the broad-host-range cloning vector pBBR1MCS, carrying different antibiotic-resistance cassettes. *Gene* **166**: 175-176.
- Tauch A, Kirchner O, Löffler B, Gotker S, Puhler A, Kalinowski J (2002). Efficient electrotransformation of *Corynebacterium diphtheriae* with a mini-replicon derived from the *Corynebacterium glutamicum* plasmid pGA1. *Curr Microbiol* **45**: 362-367.
- Xin Y, Liu H, Cui F, Liu H, Xun L (2016). Recombinant Escherichia coli with sulfide: Quinone oxidoreductase and persulfide dioxygenase rapidly oxidizes sulfide to sulfite and thiosulfate via a new pathway. *Environmental microbiology*.
- Zhang PY, Xu PP, Xia ZJ, Wang J, Xiong J, Li YZ (2013). Combined treatment with the antibiotics kanamycin and streptomycin promotes the conjugation of *Escherichia coli*. *FEMS microbiology letters* **348**: 149-156.