


Supplementary Data



SUPPLEMENTARY FIG. S1. Morphometric data at 2 dpi and 2 and 5 wpi with treatment of miR-239B (control), miR-29B delivered using a hyaluronan-based hydrogel, the hyaluronan-based hydrogel alone (in the case of the hyaluronan-based hydrogel alone, measurements were performed at 2 days and 2 and 4 wpi), and saline only injection (measurements were performed at 1 and 4 wpi). Asterisk (*) indicates a statistically significant difference in EF between the groups indicated, at the time points indicated ($p < 0.05$). dpi, days postinfarction; EF, ejection fraction; miR-29B, microRNA-29B; wpi, weeks postinfarction.