

Supporting Information

Marshall et al. 10.1073/pnas.1722493115

Table S1. Detailed CMS procedure

CMS start	Day 1	Day 2	Day 3	Day 4	Day 5	Day 6	Day 7
Week 1							
8 am	Confinement for 60 min in a 4 × 5 cm tube	Confinement for 60 min in a 4 × 5 cm tube	Food restriction (50 mg of food pellets ×3 h)	Confinement for 60 min in a 4 × 5 cm tube	Confinement for 60 min in a 4 × 5 cm tube	Access to an empty bottle (2 h)	Inversion of the light/dark cycle
1 pm	Forced bath in water at 32 °C for 30 min	Confinement for 60 min in a 4 × 5 cm tube	Confinement for 60 min in a 4 × 5 cm tube	Paired housing in damp sawdust for 10 min	Total darkness for 2 h	Housing in mildly damp sawdust for 24 h	
5 pm		Water and food deprivation for 15 h			Water deprivation for 15 h		
Week 2							
8 am	Confinement for 60 min in a 4 × 5 cm tube	Confinement for 60 min in a 4 × 5 cm tube	Confinement for 60 min in a 4 × 5 cm tube	Paired housing with food restriction for 10 min	Inversion of the light/dark cycle	Inversion of the light/dark cycle	Inversion of the light/dark cycle
1 pm	Forced bath in water at 32 °C for 30 min	Paired housing in damp sawdust for 18 h	Confinement for 60 min in a 4 × 5 cm tube	Inversion of the light/dark cycle			Housing in mildly damp sawdust for 19 h
5 pm			Water and food deprivation for 15 h				
Week 3							
8 am	Confinement for 60 min in a 4 × 5 cm tube	Food restriction (50 mg of food pellets ×3 h)	Paired housing in damp sawdust (10 min)	Forced bath in water at 32 °C (30 min)	Access to an empty bottle (2 h)		Inversion of the light/dark cycle
1 pm	Confinement for 60 min in a 4 × 5 cm tube	Confinement for 60 min in a 4 × 5 cm tube	Confinement for 60 min in a 4 × 5 cm tube	Confinement for 60 min in a 4 × 5 cm tube	Confinement for 60 min in a 4 × 5 cm tube	Inversion of the light/dark cycle	
5 pm	Water and food deprivation for 15 h			Water deprivation for 15 h	Housing in mildly damp sawdust for 20 h		