


Supplementary Information

A comprehensive overview of *FCGR3A* gene variability by full-length gene sequencing including the identification of V158F polymorphism

N.M. Mahaweni^{1,2}, T.I. Olieslagers¹, I. Olivares Rivas¹, S. J. J. Molenbroeck¹, M. Groeneweg¹, G. M. J. Bos², M. G. J. Tilanus¹, C. E. M. Voorter¹, L. Wieten¹

¹ Department of Transplantation Immunology, Tissue Typing Laboratory, GROW School for Oncology and Developmental Biology, Maastricht University Medical Center+, the Netherlands

² Department of Internal Medicine, division of Hematology, GROW School for Oncology and Developmental Biology, Maastricht University Medical Center+, the Netherlands



Supplementary Figure S1. The detection of V/F polymorphism by SSP showing 3 different genotypes.

Three different DNA samples representing homozygous T (a), heterozygous (b), and homozygous G (c) genotype. T on the lane means the SSP used for amplification was specific for the T allele, while a G indicates the SSP used for amplification was specific for the G allele. Each figure was derived from 3 different gels. Red dotted lines marked the border between the ladder and the samples.