



1 Active TSS	TssA
2 Flanking active TSS	TssAFlnk
3 Transcr. at gene 5' and 3'	TxFlnk
4 Strong transcription	Tx
5 Weak transcription	TxWk
6 Genic enhancers	EnhG
7 Enhancers	Enh
8 ZNF genes + repeats	ZNF/Rpts
9 Heterochromatin	Het
10 Bivalent/poised TSS	TssBiv
11 Flanking bivalent TSS/Enh	BivFlnk
12 Bivalent enhancer	EnhBiv
13 Repressed Polycomb	ReprPC
14 Weak repressed Polycomb	ReprPCWk
15 Quiescent/low	Quies