

Supplemental Table 2: Limits of detection for each analyte investigated from Human Angiogenesis High Performance Assay (R&D Systems, Minneapolis, MN)

Analyte	Pro or Anti-Angiogenic	Low standard of detection (pg/mL)	High standard of detection (pg/mL)
Angiogenin	Pro	2.5	10200
Angiopoietin-1	Pro	33	135,000
Endostatin	Anti	88	90,000
FGF acidic	Pro	5	20,600
FGF basic	Pro	24	5,750
PDGF-AA	Pro	1.8	7,400
PDGF-BB	Pro	5.1	20,700
PIGF	Pro	2.4	9,800
Thrombospondin	Anti	204	209,000
VEGF	Pro	4	2850
VEGF-D	Pro	18	75,000