

Run	Sample type	BMP Volatile Solids Ratio (Digestate VS: EFB VS)	pH	Feedstock added (g)	5M acetic acid added (ml)
1	Microwave EFB	2:1	7	64.35	-
1	Pressure cooked EFB	2:1	7	7.69	-
1 3	Steamed EFB	2:1	7	8.13 5.66	-
1	Microwave EFB	2:1	6	64.35	4.2
1	Pressure cooked EFB	2:1	6	7.69	5
1	Steamed EFB	2:1	6	8.13	5.1
1	Microwave EFB	1:2	7	257.4	-
1 3	Steamed EFB	1:2	7	32.51 22.64	-
3	EFB	1:2	7	12.59	-
2	Composted EFB	2:1	7	5.18	-
2	Ionic liquid EFB	2:1	7	28.89	-
2 3	EFB	2:1	7	3.15 3.15	-
2 3	Milled EFB	2:1	7	2.98 3.03	-
2	Composted EFB	2:1	6	5.18	4
2	Ionic liquid EFB	2:1	6	28.89	4
2	Untreated EFB	2:1	6	3.15	4
2 3	Digestate only	-	7	-	-