

ESM Table 1 (Antibody)

Antibody	Host	Dilution for WB	Dilution for IHC	Supplier (Antibody validated by the supplier unless indicated)
hnRNP F	Rabbit polyclonal	1:10000	1:200	specifically recognizing hnRNP F (CTARRYIGIVKQAGLER) were generated and validated in our labo ²⁵
ACE-2	Goat polyclonal	1 :2000	1 :100	R&D Systems
MasR (MAS1)	Rabbit polyclonal	1 :2000	1 :100	Novus biologicals
ACE	Goat polyclonal	1:2000	1:100	Santa Cruz Biotechnology
TGFβ1	Rabbit polyclonal	1 :2000	1 :100	Santa Cruz Biotechnology
TGF beta Receptor II	Rabbit polyclonal	1:2000	1:100	Santa Cruz Biotechnology
TGF beta Receptor I	Rabbit polyclonal	1 :2000	1 :100	Novus biologicals
Collagen type IV	Mouse Monoclonal		1 :100	Millipore
Fibronectin I	Rabbit polyclonal	1 :2000	1 :100	Sigma-Aldrich
Collagen type I	Rabbit polyclonal		1 :100	Abcam
Agt	Rabbit polyclonal	1:2000	1:200	specifically recognizing Agt were generated and validated in our lab [1]
p-Smad 2/3	Rabbit polyclonal	1 :2000		Cell signaling
Smad 2/3	Rabbit polyclonal	1 :2000		Chemicon
β-Actin	Mouse Monoclonal	1 :20000		Sigma-Aldrich

1. Wang L, Lei C, Zhang SL et al. Synergistic effect of dexamethasone and isoproterenol on the expression of angiotensinogen in immortalized rat proximal tubular cells. *Kidney international* 1998; 53: 287-295