


S2: EnzyChrom Triglyceride Assay



Diet	2020 hypoprotein	2018 cricket	milk	peanut
mmol/l	1.2460	0.5522	0.7224	0.7615
	1.6085	2.5494	1.7278	0.4889
	1.2075	0.8854	0.6767	1.1651
	0.6350	1.8002	0.3676	1.2525
	0.8378	1.0901	0.6004	0.7986
	0.6030	2.0317	1.0758	1.9906
	0.6135	1.2381	0.7778	1.1957
	1.0686	0.7973	0.5059	0.8886
	1.0947	1.5081	0.5743	1.7774
	1.2147	1.9085	0.4993	1.4018
			0.8006	0.9890
			0.7184	