Supplementary material BMJ Open

## **Supplementary Information**





Figure S1. Illustration of the applied loads via mask straps.

1 2

Supplementary material BMJ Open





27 28


29

30

**Figure S2.** Mechanical testing on samples cut directly from masks exposed to a variety of sterilization methods including 10 cycles of autoclaving, 10-minute soak in 10% bleach solution, and 10-minute soak in isopropanol.

Page 2 of 5

Supplementary material BMJ Open



**Figure S3.** Front view of 3D facial reconstruction of participants faces in fit trial of the iMASC system.

31

32 33

34

Page 3 of 5

**Figure S4.** Side view of 3D facial reconstruction of participants faces in fit trial of the iMASC system.

46

Table S1. Array of N95 and N95-comparable technologies.

Table S1. Array of N95 and Type	Ex SUS Dollars (USD)/Unit	Pros	Cons	Recommended Sterilization Method
Disposable FFR	3M 8210Plus \$1.02 USD mcmaster.com/5450T42	<ul> <li>Ease of fit/use</li> <li>Low cost per use compared to HFRR and FFRR</li> <li>Some models come with exhaust valve</li> </ul>	<ul> <li>Not reusable</li> <li>If exhaust valve is available, it's not filtered</li> <li>No eye protection</li> </ul>	N/A
iMASC system	Mask: < \$7.00 USD Filter: \$0.50 USD	<ul> <li>Low cost</li> <li>Ease/accessibility of manufacturing</li> <li>Potentially autoclavable</li> </ul>	<ul> <li>No exhaust valve for humidity/ease of use relief</li> <li>No eye protection</li> </ul>	Autoclave, Clorox wipe, IPA wipe, detergent and sterilization agent wash
Half-Face Reusable (HFRR)	Mask: \$58.98 USD mcmaster.com/5541T16- 5541T162 Replacement Cartridge: \$14.32/pair USD mcmaster.com/54445T229 Replacement Filter: \$2.65/pair USD mcmaster.com/54445T189	<ul> <li>Powered air compatible with select models</li> <li>Exhaust valve reduces humidity and breathing resistance</li> <li>Flange/gusset provides comfortable seal skin</li> </ul>	<ul> <li>Expensive</li> <li>Exhaust valve not filtered</li> <li>No eye protection</li> </ul>	Detergent and sterilization agent wash
Full-Face Reusable (FFRR)	Mask: \$168.47 USD mcmaster.com/5541T28 Replacement Filter: \$2.65/pair USD mcmaster.com/54445T189	<ul> <li>Best coverage</li> <li>Powered air compatible with select models</li> <li>Exhaust valve reduces humidity and breathing resistance</li> <li>Comfortable seal to skin</li> </ul>	<ul> <li>Expensive</li> <li>Exhaust valve not filtered</li> <li>Potential visual obstruction due to fogging</li> </ul>	Detergent and sterilization agent wash