

Supplementary Table 2. List of plasmids used for the study.

Gene	Backbone	Tag	Mutation	Resource	Reference
CFL1	pCMV6-Entry	Flag-Myc	-	ORIGENE	-
CHCHD3	pCMV6-Entry	Flag-Myc	-	ORIGENE	-
EFHD2	pCMV6-Entry	Flag-Myc	-	ORIGENE	-
KHDRBS1	pCMV6-Entry	Flag-Myc	-	ORIGENE	-
LASP1	pCMV6-Entry	Flag-Myc	-	ORIGENE	-
LIMA1	pCMV6-Entry	Flag-Myc	-	ORIGENE	-
PPIL3	pCMV6-Entry	Flag-Myc	-	ORIGENE	-
STEEP	pCMV6-AC	-	-	ORIGENE	-
TBK1	pcDNA3	Flag	-	Soren R Paludan	-
TRIF	pcDNA3	-	-	Soren R Paludan	-
MAVS	pcDNA3	-	-	Soren R Paludan	-
cGAS	pcDNA3	Flag	-	Soren R Paludan	-
IRF3	pcDNA3		S396D, S398D, S402D, T404D, and S405D	Soren R Paludan	-
STING	pCMV-Tag 4	Flag	-	Soren R Paludan	-
STEEP	pCMV6-AN-HA	HA	-	Soren R Paludan	-
STEEP	pCMV6-AN-HA	HA	R7A	Soren R Paludan	-
STEEP	pCMV6-AN-HA	HA	Y30A,C34A	Soren R Paludan	-
STEEP	pCMV6-AN-HA	HA	R51A,D54A,V58A	Soren R Paludan	-
STEEP	pCMV6-AN-HA	HA	E86A, C93A,C96A	Soren R Paludan	-
STEEP	pCMV6-AN-HA	HA	K149A,G152A,E166A, E168A,Y178A	Soren R Paludan	-
Ubiquitin	pRK5-HA	HA	K11 only	Ted Dawson	1
Ubiquitin	pRK5-HA	HA	K27 only	Ted Dawson	1
Ubiquitin	pRK5-HA	HA	K29 only	Ted Dawson	1
Ubiquitin	pRK5-HA	HA	K33 only	Ted Dawson	1
Ubiquitin	pRK5-HA	HA	K48 only	Ted Dawson	1
Ubiquitin	pRK5-HA	HA	K63 only	Ted Dawson	1
hVps34	pVitro2	His, V5	-	Jonathan Backer	2
hVps15					
pCI-neo-HA-hAtg14	pCI-neo HA Ci3	HA		Noboru Mizushima	3
UVRAG	pRK5	HA		Do-Hyung Kim	4
Beclin 1	pRK5	HA		Do-Hyung Kim	5
Beclin 1	pCDNA4 - 3xflag	Flag		Qing Zhong	6
SAR1	pCMV6-Entry	Flag	-	ORIGENE	-
SAR1	pCMV6-Entry	Flag	H79G	Soren R Paludan	-
GFP ¹³³	pEGFP-N1	GFP	-	Christine M. Doucet	7
GFP-FYVE	peGFP-C2	GFP	-	Harald Stenmark	8
MTMR3	mStrawberry-C2	mCherry	-	Takeshi Noda	9
MTMR3	mStrawberry-C2	mCherry	C143S	Takeshi Noda	9
STING	pcDNA3	Flag	-	Zhijian J Chen	10
STING	pcDNA3	Flag	R238A,240A	Zhijian J Chen	10
STING	pcDNA3	Flag	Q273A,Q277A	Zhijian J Chen	10
STING	pcDNA3	Flag	L374A, R375A	Zhijian J Chen	11
STING	pcDNA3	Flag	V155M	Zhijian J Chen	11
STING	pcDNA3	Flag	G158A	Zhijian J Chen	11
VSVG	pcDNA3	mCherry	-	Jacob G Mikkelsen	-

ATG14	pCMV3-GFPSpark	GFP	-	Sino Biological	-
UVRAG	pCMV3-GFPSpark	GFP	-	Sino Biological	-

1. **Parkin mediates nonclassical, proteasomal-independent ubiquitination of synphilin-1: implications for Lewy body formation.** Lim KL, Chew KC, Tan JM, Wang C, Chung KK, Zhang Y, Tanaka Y, Smith W, Engelender S, Ross CA, Dawson VL, Dawson TM. *J Neurosci*. 2005 Feb 23; 25(8):2002-9. 10.1523/JNEUROSCI.4474-04.2005 [PubMed 15728840](#)
2. **hVps15, but not Ca2+/CaM, is required for the activity and regulation of hVps34 in mammalian cells.** Yan Y, Flinn RJ, Wu H, Schnur RS, Backer JM. *Biochem J*. 2009 Feb 1. 417(3):747-55. 10.1042/BJ20081865 [PubMed 18957027](#)
3. **mTORC1 phosphorylates UVRAG to negatively regulate autophagosome and endosome maturation.** Kim YM, Jung CH, Seo M, Kim EK, Park JM, Bae SS, Kim DH. *Mol Cell*. 2015 Jan 22;57(2):207-18. doi: 10.1016/j.molcel.2014.11.013. Epub 2014 Dec 18. 10.1016/j.molcel.2014.11.013 [Pubmed 25533187](#)
4. **Beclin 1 forms two distinct phosphatidylinositol 3-kinase complexes with mammalian Atg14 and UVRAG.** Itakura E, Kishi C, Inoue K, Mizushima N. *Mol Biol Cell*. 2008 Dec . 19(12):5360-72. 10.1091/mbc.e08-01-0080 [PubMed 18843052](#)
5. **The ULK1 complex mediates MTORC1 signaling to the autophagy initiation machinery via binding and phosphorylating ATG14.** Park JM, Jung CH, Seo M, Otto NM, Grunwald D, Kim KH, Moriarity B, Kim YM, Starker C, Nho RS, Voytas D, Kim DH. *Autophagy*. 2016;12(3):547-64. doi: 10.1080/15548627.2016.1140293. 10.1080/15548627.2016.1140293 [PubMed 27046250](#)
6. **Identification of Barkor as a mammalian autophagy-specific factor for Beclin 1 and class III phosphatidylinositol 3-kinase.** Sun Q, Fan W, Chen K, Ding X, Chen S, Zhong Q. *Proc Natl Acad Sci U S A*. 2008 Dec 9. 105(49):19211-6. 10.1073/pnas.0810452105 [PubMed 19050071](#)
7. Doucet, Christine M., et al. "Membrane curvature sensing by amphipathic helices is modulated by the surrounding protein backbone." *PloS one* 10.9 (2015): e0137965.
8. Raiborg, Camilla, et al. "Repeated ER–endosome contacts promote endosome translocation and neurite outgrowth." *Nature* 520.7546 (2015): 234.
9. Taguchi-Atarashi, Naoko, et al. "Modulation of local PtdIns3P levels by the PI phosphatase MTMR3 regulates constitutive autophagy." *Traffic* 11.4 (2010): 468-478.
10. Shang, Guijun, et al. "Cryo-EM structures of STING reveal its mechanism of activation by cyclic GMP–AMP." *Nature* 567.7748 (2019): 389.
11. Zhang, Conggang, et al. "Structural basis of STING binding with and phosphorylation by TBK1." *Nature* 567.7748 (2019): 394.