

Biohybrid Bovine Bone Matrix for Controlled Release of Mesenchymal Stem/Stromal Cell Lyosecretome: a Device for Bone Regeneration

Elia Bari ^{1#}, Ilaria Roato ^{2#}, Giuseppe Perale ^{3,4,5}, Filippo Rossi ⁶, Tullio Genova ⁷, Federico Mussano ², Riccardo Ferracini ⁸, Marzio Sorlini ^{9,10}, Maria Luisa Torre ^{1,10*}, Sara Perteghella ^{1,10}

¹ Department of Drug Sciences, University of Pavia, Viale Taramelli 12, I-27100 Pavia, Italy; elia.bari@unipv.it (E.B.); sara.perteghella@unipv.it (S.P.); marina.torre@unipv.it (M.L.T.)

² CIR-Dental School, Department of Surgical Sciences, University of Torino, Via Nizza 230, 10126 Torino Italy; ilaria.roato@unito.it (I.R.); federico.mussano@unito.it (F.M.)

³ Industrie Biomediche Insubri SA, via Cantonale 67, 6805 Mezzovico-Vira, Switzerland; giuseppe@ibi-sa.com (G.P.)

⁴ Faculty of Biomedical Sciences, University of Southern Switzerland (USI), Via G. Buffi 13, 6900 Lugano, Switzerland

⁵ Ludwig Boltzmann Institute for Experimental and Clinical Traumatology, Donaueschingenstrasse 13, 1200 Vienna, Austria

⁶ Department of Chemistry, Materials and Chemical Engineering "Giulio Natta", Politecnico di Milano, via Mancinelli 7, 20131 Milano, Italy; filippo.rossi@polimi.it

⁷ Department of Life Sciences and Systems Biology, University of Torino, Via Accademia Albertina 13, 10123 Torino, Italy; tullio.genova@unito.it

⁸ Department of Surgical Sciences and Integrated Diagnostics, University of Genova, Viale Benedetto XV 6, 16132 Genova, Italy; ferracini@edu.unige.it

⁹ SUPSI - Department of Innovative Technologies, Lugano University Centre, Campus Est, Via la Santa 1, 6962 Viganello, Switzerland; marzio.sorlini@supsi.ch

¹⁰ PharmaExceed Srl, Piazza Castello 19, I-27100 Pavia, Italy

* Correspondence: marina.torre@unipv.it; Tel.: +390382987779

Equally contributed to the work.




Figure S1. Morphological investigation by ESEM of SmartBone[®] loaded with 0.1% *w/v* poloxamer 407 (a-c). SEM images were taken at increasing magnifications (30× and 100×). Scale bar: 500 and 100 µm.