

A

MiSeq sequencing and data analysis of CD46-Cas9 cells and controls


Outknocker web tool

CD46 Cas9 cells:

Knockout efficiency: 84.2 %

WT cells:

Knockout efficiency: 0.1 %




REFERENCE CTATGGAGCTCATTGGTAAACCAAAAC**CCT**ACTATGAGATTGGTGAACGAGTAGATTATAAGTGTAAAAAAGGATACTTC
PAM

Intact		15.71%
no indel	CTATGGAGCTCATTGGTAAACCAAAAC CCT ACTATGAGATTGGTGAACGAGTAGATTATAAGTGTAAAAAAGGATACTTC	919 reads
Deleted		74.35%
90nt deletion	CTATGGAGCTCATTGGTAAACCAAAAC CCT TACTAT	274 reads
91nt deletion	CTATGGAGCTCATTGGTAAACCAAAAC CCT TACT	2995 reads
92nt deletion	CTATGGAGCTCATTGGTAAACCAAAAC CCT TACTA	1081 reads
Indel		9.94%
1nt deletion	CTATGGAGCTCATTGGTAAACCAAAAC CCT TACT_TGAGATTGGTGAACGAGTAGATTATAAGTGTAAAAAAGGATACTTC	316 reads
2nt deletion	CTATGGAGCTCATTGGTAAACCAAAAC CCT TACT__GAGATTGGTGAACGAGTAGATTATAAGTGTAAAAAAGGATACTTC	265 reads

B

non-infected



EBV-infected

