

Supplement Table 2. Semen characteristics of patients with follow-up

Variable	Patient No.				
	1	2	3	4	5
Semen characteristics following COVID-19 diagnosis					
Volume (mL)	2.6	2.2	2.0	2.0	2.1
Sperm concentration ($\times 10^6$ /mL)	22	20	11	5	0.17
Total sperm number (million)	57.2	44	22	10	0.36
Sperm total motility (%)	0	0	0	0	14
Semen (pH)	7.6	7.4	7.6	6.8	7.6
Follow-up semen characteristics					
Volume (mL)	0.3	1.5	1.6	2.1	1.0
Sperm concentration ($\times 10^6$ /mL)	20	12	21	20	12
Total sperm number (million)	6	18	33.6	42	12
Sperm total motility (%)	0	0	0	69	39
Semen (pH)	7.6	7.6	7.2	7.6	7.6
In-person/mailed sample ^a	Mail	Mail	Mail	In-person	In-person

Values are presented as median.

COVID-19: coronavirus disease.

^aIn-person: samples provided in-office, mail: samples were shipped overnight.