

| | | 1 cycle 40 s | | | 2 cycles of 40 s | | | 3 cycles of 40 s | | | 4 cycles of 30 s | | |
|--------------------------------------|--------------------|--------------|-------|-------|------------------|-------|-------|------------------|-------|-------|------------------|-------|-------|
| <i>Strongyloides stercoralis</i> | FastPrep-24™ 5G | 27.72 | 27.59 | 27.53 | 26.9 | 26.66 | 26.99 | 27.19 | 27.16 | 27.3 | 26.08 | 25.87 | 26.11 |
| | Mini-Beadbeater-24 | 27.05 | 26.96 | 26.96 | 27.67 | 27.66 | 27.53 | 26.62 | 27.49 | 27.74 | 26.85 | 26.78 | 26.98 |
| | | | | | | | | | | | | | |
| <i>Ascaris lumbricoides</i> | FastPrep-24™ 5G | 32.62 | 32.66 | 32.97 | 32.42 | 32.1 | 31.43 | 32.54 | 32.01 | 31.68 | 32.31 | 32.26 | 32.56 |
| | Mini-Beadbeater-24 | 32.18 | 32.34 | 32.35 | 32.69 | 32.9 | 32 | 33.7 | 32.9 | 33.26 | 32.36 | 32.88 | 32.76 |