## **Biomedical Optics EXPRESS**

## Spatially resolved measurements of ballistic and total transmission in microscale tissue samples from 450 nm to 1624 nm: supplement

AARON T. MOK,<sup>1,2,\*</sup> D JAMIEN SHEA,<sup>3</sup> CHUNYAN WU,<sup>1,4</sup> D FEI XIA,<sup>1,2,5</sup> D ROSE TATARSKY,<sup>3</sup> NILAY YAPICI,<sup>3</sup> AND CHRIS XU<sup>1</sup>

This supplement published with Optica Publishing Group on 21 December 2021 by The Authors under the terms of the Creative Commons Attribution 4.0 License in the format provided by the authors and unedited. Further distribution of this work must maintain attribution to the author(s) and the published article's title, journal citation, and DOI.

Supplement DOI: https://doi.org/10.6084/m9.figshare.17147687

Parent Article DOI: https://doi.org/10.1364/BOE.441844

<sup>&</sup>lt;sup>1</sup>School of Applied and Engineering Physics, Cornell University, NY 14853-3501, USA

<sup>&</sup>lt;sup>2</sup>Meining School of Biomedical Engineering, Cornell University, NY 14853-3501, USA

<sup>&</sup>lt;sup>3</sup>Department of Neurobiology and Behavior, Cornell University, NY 14853-3501, USA

<sup>&</sup>lt;sup>4</sup>College of Veterinary Medicine, Cornell University, NY 14853-3501, USA

<sup>&</sup>lt;sup>5</sup>Present address: Laboratoire Kastler Brossel, ENS-PSL Research University, CNRS, Sorbonne Université, Collège de France, Paris, France

<sup>\*</sup>tm624@cornell.edu

## Spatially resolved measurements of optical properties in microscale biological samples from 450nm to 1624nm: supplemental document




Fig. S1. Transmission measurement of  $100\mbox{-}\mu m$  thick UV curable resin using a commercial spectrophotometer.