

Supporting Information

A Highly Efficient Nano-Sized Cu₂O/SiO₂ Egg-Shell Catalyst for C-C Coupling Reactions

Soohee Kim,^a Shin Wook Kang,^b Ji Chan Park,^{*b} Kang Hyun Park^{*a}

^a Department of Chemistry and Chemistry Institute for Functional Materials, Pusan National University, Busan 609-735, Korea.

Phone: (+82)-51-510-2238; Fax: (+82)-51-980-5200; e-mail: chemistry@pusan.ac.kr

^b Clean Fuel Laboratory, Korea Institute of Energy Research, Daejeon, 305-343, Korea. e-mail: jcpark@kier.re.kr

Entry	Solvent (2 ml)	Conv. (%)	Select. (%)
1	Toluene	84	99
2	THF	3	100
3	DMF	19	0
4	neat	100	99

Reaction conditions: phenylacetylene (0.50 mmol), benzoyl chloride (0.75 mmol), base Et₃N (4.0 equiv.), cat. 1 mol% SiO₂@Cu₂O/mSiO₂, 12h, 80 °C, Ar atmosphere. Determined by using gas chromatography-mass spectrometry (GC-MS).

Table S1 The synthesis of ynones from acyl chlorides and terminal alkynes catalyzed by SiO₂@Cu₂O/mSiO₂ under the various solvent system.