



**Figure S5. Related to Fig. 7. Inhibition of clot retraction by  $\alpha\text{IIb}\beta_3$  inhibitors.** (A) Human PRP was pre-incubated with buffer control or indicated concentrations of  $\alpha\text{IIb}\beta_3$  inhibitors, and then thrombin was added at 0.75 U/ml to stimulate clot retraction. The tubes were imaged at 30 and 90 mins. (B) Clot retraction was quantified with ImageJ and presented as 1- (clot area/whole reaction area). Data were fitted with the function of inhibitor vs response (three parameters) in Prism 8.