

Supporting Information

Enhancement of Magnetic Surface-Enhanced Raman Scattering Detection by Tailoring Fe₃O₄@Au Nanorod Shell Thickness and its Application in the On-site Detection of Antibiotics in Water

Leixuri B. Berganza,¹ Lucio Litti,^{2} Moreno Meneghetti,² Senentxu Lanceros-Méndez,^{1,3}*

Javier Reguera,^{1}*

¹BCMaterials, Basque Center for Materials, Applications, and Nanostructures, UPV/EHU
Science Park, Leioa 48940, Spain.

²Nanostructures and Optics Laboratory, Department of Chemical Sciences, University of
Padova, Via Marzolo, 1, 35131 Padova, Italy.

³Ikerbasque Basque Foundation for Science Bilbao 48009, Spain.

J. Reguera: javier.reguera@bcmaterials.net; *L. Litti:* lucio.litti@unipd.it

KEYWORDS: Nanorod, plasmonic, magnetic, SERS, antibiotics, portable Raman,
Ciprofloxacin




Figure S1. Low magnification TEM image of Fe₃O₄-NR: Au core-satellite nanoparticles used as seed for further Au growth.




Figure S2. Low magnification TEM image of Au coated nanorods corresponding to sample NR@Au1.




Figure S3. Low magnification TEM image of Au coated nanorods corresponding to sample NR@Au2.




Figure S4. Low magnification TEM image of Au coated nanorods corresponding to sample NR@Au3.




Figure S5. Low magnification TEM image of Au coated nanorods corresponding to sample NR@Au4.




Figure S6. Low magnification TEM image of Au coated nanorods corresponding to sample NR@Au5.




Figure S7. Low magnification TEM image of Au coated nanorods corresponding to sample NR@Au6