

Additional file 1

List of primary and secondary antibodies.

Primary antibodies	Supplier	Species raised in; clonality	Dilution (application)	Reference
CXCR4	Abcam, Cambridge, UK	Goat, polyclonal	1:1000 (WB)	ab1670
CXCR7	Acris, Herford, Germany	Rabbit, polyclonal	1:1000 (WB)	AP17961PU-N
CXCR3	R&D Systems, Minneapolis, MN	Mouse, monoclonal	1:1000 (WB)	clone #49801
cleaved caspase-3 (N175)	Cell Signaling Technologies, Danvers, MA	Rabbit, polyclonal	1:400 (ICC)	#9661
GAPDH	Fitzgerald Industries International, Acton, MA	Mouse, monoclonal	1:5000 (WB)	10RG109A
Secondary antibodies	Supplier	Species raised in; clonality	Dilution (application)	Reference
Alexa Fluor 488 anti-rabbit	Thermo Fisher Scientific, Waltham, MA	Donkey, polyclonal	1:200 (ICC)	A21206
Alexa Fluor 555 anti-rabbit	Thermo Fisher Scientific	Donkey, polyclonal	1:200 (ICC)	A31572
Peroxidase-conjugated IgG anti-rabbit	Jackson ImmunoResearch, West Grove, PA	Donkey, polyclonal	1:10,000 (WB)	711-035-152
Peroxidase-conjugated IgG anti-mouse	Vector Laboratories, Burlingame, CA	Horse	1:10,000 (WB)	PI-2000

WB, Western blot; ICC, immunocytochemistry